Stat Profile Prime Plus® Calibrator Cartridge With Creatinine

Stat Profile Prime Plus® Kalibrerkassette mit Kreatinin, Stat Profile Prime Plus® Φυσίγγιο μέσου βαθμονόμησης με κρεατινίνη, Cartucho de calibrador con creatinina Stat Profile Prime Plus®, Cartouche de calibrateurs Stat Profile Prime Plus® avec créatinine, Cartuccia calibratore con creatinina Stat Profile Prime Plus®, Cartucho de calibrador Stat Profile Prime Plus® com creatinina, Stat Profile Prime Plus® kalibrátor patron kreatininnel, კრეატინის უსაფრთხო გამაზნი Stat Profile Prime Plus®, Stat Profile Prime Plus® キャリブレータートリッジ (クリアチニン), 캐리브레이터 카트리지, Stat Profile Prime Plus® 보정기 카트리지, Stat Profile Prime Plus® 含肌酐的校准试剂盒

LOT 23152023

2024-11-09

CAL A	CAL B	CAL C	CAL D	CAL F	SOLN R
HCO ₃	11 mmol/L	pH	6.834 @ 37°C	HCO ₃	22 mmol/L
Na ⁺	135.5 mmol/L	H ⁺	146.55 nmol/L	Cl ⁻	69.0 mmol/L
K ⁺	4.0 mmol/L	Na ⁺	77.0 mmol/L	Glu	200 mg/dL
Cl ⁻	108.0 mmol/L	K ⁺	10.0 mmol/L	Glu	11.10 mmol/L
iCa	1.04 mmol/L	iCa	1.92 mmol/L	Lac	10.0 mmol/L
iCa	4.16 mg/dL	iCa	7.68 mg/dL	Lac	89.0 mg/dL
iMg	0.50 mmol/L	iMg	1.50 mmol/L	BUN	50 mg/dL
iMg	1.22 mg/dL	iMg	3.63 mg/dL	BUN	17.9 mmol/L
Glu	80 mg/dL	tHb	15.0 g/dL	Creat	5.0 mg/dL
Glu	4.44 mmol/L	tHb	150 g/L	Creat	0.44 mmol/L
Lac	2.0 mmol/L	tHb	9.3 mmol/L	Creat	442 µmol/L
Lac	17.8 mg/dL				
BUN	10 mg/dL				
BUN	3.6 mmol/L				
Creat	1.0 mg/dL				
Creat	0.09 mmol/L				
Creat	88.4 µmol/L				

(EN)

Description

This Calibrator contains the following aqueous solutions within individual bags: 8 Calibrators (Calibrators A, B, C, D and F), Reference Solution, and a flexible waste bag. Each bag includes a filament with a septum. The exposed bag filaments are arranged in a line along the rear of the cartridge box. The septa are pierced during the insertion of the cartridge into the analyzer. A charge kit consisting of 2 syringes filled with creatinine solution is included for creatinine activation prior to installation.

Intended Use

For the quantitative determination of Na⁺, K⁺, Cl⁻, iCa, iMg, BUN (Urea), Glucose, Lactate, Creatinine, pH, PO₂, PO₃, SO₄, Hct, tHb, HbF, O₂Hb, COHb, MetHb, HHb, and tBil¹ in human blood using the Stat Profile Prime Plus Analyzer.

Methodology

Refer to Stat Profile Prime Plus Analyzer Instructions For Use Manual for Methodologies and Principles.

Composition

The Calibrator contains aqueous salt solutions containing analytes at concentrations listed in the table above. Calibrator B also contains dyes. Contain no constituents of human origin, however, good laboratory practice should be followed during the handling of these materials. (REF. NCLCS DOCUMENT M29-T2).

Warnings and Cautions

Intended for in vitro diagnostic use.

Refer to the Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory reagents. Discard Cartridge in accordance with local guidelines.

Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

Storage

Store at 15 - 30 °C. Expiration date is printed on the Calibrator Cartridge.

Directions For Use**Activation of Cartridge**

Cal A and Cal C must be charged with creatinine syringe prior to installation of the Calibrator Cartridge. Syringes are labeled and color coded to correspond to the filaments on the cartridge for Cal A and Cal C.

1. Hold the syringe with tip side down and remove protective cap.

2. Attach one of the enclosed needles assemblies to the syringe. Remove the protective cover from the needle.

3. Gently tap the cap and label of the syringes to the appropriate filament and insert needle.

4. Slowly depress syringe plunger until the contents are dispensed.

DO NOT PUSH BACK ON THE PLUNGER TO FLUSH CONTENTS OF SYRINGE.

5. Remove needle/syringe assembly from filament and discard in an appropriate sharps container.

6. Repeat Steps 1-5 for the next Cal.

7. Mix Cartridge well by gently shaking for 1 minute.

Refer to the Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete directions. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents. Cartridge is ready for use..

*Not available in the USA or for Point-of-Care/Near-Patient Testing use.

Traceability of Standards

Analyses traced to NIST Standard Reference Materials wherever possible.

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