Stat Profile® pHox® Ultra/pHox/Critical Care Xpress SO₂ Calibrator

SO₂ Kalibrator, Βαθμονομητής SO₂, Calibrador de SO₂, Etalon SO₂, Calibratore SO₂,
SO₂ キャリブレータ, Calibrador de SO₂, SO₂ Kalibrator

LOT 23157082

2024-11-17

SO₂ CALIBRATOR 1 | 2

Calibrator Value, Kalibratorwert, Τιμή βαθμονομητή, Valor del calibrador, Valeur étalon, Valore calibratore, キャリブレータ値, Valor do calibrador, Kalibratorvärde			
CALIBRATOR 1	LOT 23137031	2024-11-17	52.6
CALIBRATOR 2	LOT 23137032	2024-11-17	97.8

Product Description

An aqueous material for calibrating the oxygen saturation channel for use with Nova Biomedical analyzers ONLY.

Intended Use

For *in vitro* diagnostic use for calibrating SO₂% on the Nova Biomedical Stat Profile pHox/Plus/LC Analyzers, Stat Profile pHox Ultra Analyzer, and Stat Profile CCX Analyzers

Methodology

Refer to the Stat Profile pHox Ultra, pHox/Plus/LC, and STP CCX Analyzer Instructions For Use Manuals.

Composition

An aqueous electrolyte solution whose reflectance characteristics give a signal that is equivalent to a known oxygen saturation value in whole blood. Each ampule contains a minimum volume of 1.0mL. Contains no constituents of human origin, however, good laboratory practice should be followed during the handling of these materials. (REF: NCCLS DOCUMENT M29-T2).

Warnings and Precautions

Intended for *in vitro* diagnostic use.
Refer to the Stat Profile pHox/Plus/LC Analyzer and Stat Profile Critical Care Xpress Analyzer Instructions for Use Manual for complete directions for use, including recommendations for use of controls, troubleshooting information, and methodology and Principles of the testing procedures. Follow standard practices required for handling laboratory reagents.

Storage

Store at 15-30°C. DO NOT FREEZE. The expiration date is printed on each ampule.

Directions For Use

Contents must be shaken well before use. Refer to analyzer Instructions for Use Manual for complete directions. Verify that the lot number appearing on the Calibrator table is the same lot number indicated on the ampule.

Limitations

The Oxygen Saturation External Calibrator is intended for calibrating Nova Biomedical Analyzers ONLY. The external calibrator is not suitable for use on any other manufacturers' analyzer.

Traceability of Standards

Analyses are traceable to tonometry.

Assigned/Calibrator Values

The CALIBRATOR VALUE for each level was determined at Nova Biomedical by using multiple runs of each level of calibrator at 37°C on multiple instruments.

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