REF 48921, 48922

IVD 15°C 30°C CE

Stat Profile® pHOx Ultra/Critical Care Xpress CO-Oximeter Controls Auto-Cartridge

CO-Oximeter-Kontrollen Auto-QK-Kassette, Autómatr Kaoéta Díalumuráton Eléggou Oázumetrou-CO, Cartucho automático para controles de CO-Oxímetro, Cartouche CQ d'oxymètre de CO, Cartuccia automatica controlli qualità co-ossimetro, Cartucho automático de controlo de co-oxímetro, Autokassett für CO-oximeter-kontroller, CO-酸素濃度計コントロール自動カートリッジ, CO- 산소 농도계 조절 자동 카트리지, CO-血氧计质控品自动筒

CONTROL | 7 | 8 | 9

LOT 23159093 2025-04-11

Expected Ranges, Erwartete Bereiche, Αναμενόμενη περιοχή τιμών, Rangos esperados, Fourchettes attendues, Intervalli previsti, Gamas previstas, Förväntade områden, 期待範囲値, 예상 범위, 期望范围

		CONTROL 7	CONTROL 8	CONTROL 9
		min - \bar{x} - max	min - \bar{x} - max	min - \bar{x} - max
tHb	g/dL	6.8 - 7.5 - 8.2	16.0 - 17.0 - 18.0	18.1 - 19.6 - 21.1
tHb	g/L	68 - 75 - 82	160 - 170 - 180	181 - 196 - 211
tHb	mmol/L	4.2 - 4.7 - 5.1	9.9 - 10.5 - 11.2	11.2 - 12.2 - 13.1
O ₂ Hb	%	79.7 - 84.7 - 89.7	46.9 - 50.9 - 54.9	19.5 - 22.0 - 24.5
COHb	%	2.7 - 6.7 - 10.7	17.6 - 21.6 - 25.6	28.7 - 32.7 - 36.7
MetHb	%	4.0 - 7.0 - 10.0	16.0 - 19.0 - 22.0	26.7 - 29.7 - 32.7
HHb	%	0.0 - 4.0 - 8.0	4.5 - 8.5 - 12.5	11.4 - 15.4 - 19.4
tBil	mg/dL	4.1 - 4.5 - 4.9	8.6 - 10.6 - 12.6	14.3 - 18.3 - 22.3

Product Description

An assayed aqueous quality control material for monitoring the measurement of total hemoglobin (tHb), oxyhemoglobin (O₂Hb), carboxyhemoglobin (COHb), metahemoglobin (MetHb), deoxyhemoglobin (HHb) and total bilirubin (tBil). The controls are formulated at three levels.

Intended Use

For in vitro diagnostic use by healthcare professionals for monitoring the performance of Nova Biomedical Stat Profile pHOx Ultra/Critical Care Xpress Analyzers.

Methodology

Refer to Stat Profile pHOx Ultra/CCX Analyzer Instructions For Use Manual.

Composition

Contains dyes, viscosity enhancers, and preservative in an aqueous salt solution. Each control bag contains a minimum volume of 90 mL. Contains no components of human origin, however good laboratory practices should be followed during handling of these materials.

(REF: NCCLS DOCUMENT M29-T2).

Warnings and Cautions

Intended for in vitro diagnostic use. Refer to the Stat Profile pHOx Ultra/CCX Analyzer Instructions for Use Manual for complete instructions. Use standard practices for handling reagents. Methohemoglobin (MetHb) and total bilirubin (tBil) are traceable by using Cyanmethemoglobin method. Carboxyhemoglobin (COHb) and Oxyhemoglobin (O₂Hb) are traceable using Spectrophotometry.

Storage

Store at 15-30°C. Expiration date is printed on each cartridge.

Directions for use

The Auto-Cartridge should be mixed by gently inverting the pack for several seconds. DO NOT SHAKE THE PACK. Refer to Analysis Instructions for Use Manual for complete instructions. Verify that the lot number appearing on the Expected Ranges Table is the same lot number indicated on the cartridge.

Limitations

The Expected Range values are specific for instruments and calibrators manufactured by Nova Biomedical.

Traceability of Standards

Total Hemoglobin (tHb) and MetHemoglobin (MetHb) are traceable by using Cyanmethemoglobin method. Carboxyhemoglobin (COHb) and Oxyhemoglobin (O₂Hb) are traceable using Spectrophotometry.

Reference Intervals

Concentrations are formulated to represent three levels.

Users may wish to determine MEAN VALUES and EXPECTED RANGES in their own laboratory.

Expected Range

The EXPECTED RANGE for each parameter was determined at Nova Biomedical by using multiple runs of each level of control on properly calibrated instruments. The EXPECTED RANGE indicates the maximum deviations from the mean value that may be expected under differing laboratory conditions for instruments operating within specifications. Refer to Expected Ranges Table.

'How to Define and Determine Reference Intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13'

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