REF 58404, 58405

15°C 30°C CE

Stat Profile Prime Plus® VET Calibrator Cartridge With Creatinine

Stat Profile Prime Plus® VET Kalibrierkassette mit Kreatinin, Φυσιγγιο μέσου βαθμονόμησης VET με κρεατινίνη Stat Profile Prime Plus®, Cartucho de calibrador con creatinina Stat Profile Prime Plus® VET, Cartouche de calibrateurs avec créatinine Stat Profile Prime Plus® VET, Cartuccia calibratore con creatinina Stat Profile Prime Plus® VET, Cartucho de calibrador Stat Profile Prime Plus® VET com creatinina, Stat Profile Prime Plus® VET kalibrátor patron kreatininnel, Stat Profile Prime Plus® VET קריינטינין כוֹל עַם קָרִיאטִינִין כְּיֻלְּ מִיחָשְׁנִין, Stat Profile Prime Plus® VET キアリブレータ(クリアチニン), 크레이티닌 사용 Stat Profile Prime Plus® VET 보정기 카트리지, Stat Profile Prime Plus® 含肌酐的校准试剂盒

LOT 23187038

2024-12-15

CAL A	≥ 495 mL	CAL B	≥ 250 mL	CAL C	≥ 245 mL	CAL D	≥ 250 mL	CAL F	≥ 850 mL	SOLN R	≥ 500 mL
HCO ₃	11 mmol/L	pH	6.834 @ 37°C	HCO ₃	22 mmol/L	pH	7.354 @ 37°C	PO ₂	101.7 mmHg	KCl	2 mol/L
Na ⁺	135.5 mmol/L	H ⁺	146.55 nmol/L	Cl ⁻	69.0 mmol/L	H ⁺	44.26 nmol/L	PO ₂	13.53 kPa		WASTE CONTAINER W
K ⁺	4.0 mmol/L	Na ⁺	77.0 mmol/L	Glu	200 mg/dL						
Cl ⁻	108.0 mmol/L	K ⁺	10.0 mmol/L	Glu	11.10 mmol/L						
iCa	1.04 mmol/L	iCa	1.92 mmol/L	Lac	10.0 mmol/L						
iCa	4.16 mg/dL	iCa	7.68 mg/dL	Lac	89.0 mg/dL						
iMg	0.50 mmol/L	iMg	1.50 mmol/L	BUN	50 mg/dL						
iMg	1.22 mg/dL	iMg	3.63 mg/dL	BUN	17.9 mmol/L						
Glu	80 mg/dL	tHb	15.4 g/dL	Creat	5.0 mg/dL						
Glu	4.44 mmol/L	tHb	154 g/L	Creat	0.44 mmol/L						
Lac	2.0 mmol/L	tHb	9.6 mmol/L	Creat	442 µmol/L						
Lac	17.8 mg/dL										
BUN	10 mg/dL										
BUN	3.6 mmol/L										
Creat	1.0 mg/dL										
Creat	0.09 mmol/L										
Creat	88.4 µmol/L										

EN

Description

The Cartridge contains the following aqueous solutions within individual bags: 5 Calibrators (Calibrator A, B, C, D, and F), Reference Solution, and a flexible waste bag. Each bag includes a unique color-coded label. The exposed bag filaments are arranged in a line along the rear of the cartridge box. The syringes are pierced during the insertion of the cartridge into the analyzer. A charge kit consisting of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation.

Intended Use

For the quantitative determination of Na⁺, K⁺, Cl⁻, iCa, iMg, BUN (Urea), Glucose, Lactate, Creatinine, pH, PO₂, PO₃%, Hct/Hb, O₂Hb, COHb, MetHb, HHb, and iBil in animal whole blood using the Stat Profile Prime Plus VET Analyzer.

Methodology

Refer to Stat Profile Prime Plus VET Analyzer Instructions For Use Manual for Methodologies and Principles.

Composition

Calibrator and solution volumes are listed in the table. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF. NCCLS DOCUMENT M29-T2).

Warnings and Cautions

Refer to the Stat Profile Prime Plus VET Analyzer Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory reagents. Discard Cartridge in accordance with local guidelines.

Vial Bag Contains Preservative (eye, skin, and respiratory irritant).

Storage

Store at 15 - 30°C. Expiration date is printed on the Calibrator Cartridge.

Directions For Use**Activation of Cartridge**

Cal A and C must be charged with creatinine syringe prior to installation of the Calibrator Cartridge. The two syringes are labeled and color coded to correspond to the filaments on the carton for Cal A and Cal C.

- Hold the syringe with tip side down and remove protective cap.
- Connect the needleless needle assemblies to the syringes. Remove the protective caps from the needles.
- Match the color and label of the syringe to the appropriate filament and insert needle.
- Slowly depress syringe plunger until the contents are dispensed.
- DO NOT SQUEEZE BACK ON THE PLUNGER TO FLUSH CONTENTS OF SYRINGE.
- Discard the needleless assembly from filament and discard in an appropriate sharps container.
- Repeat Steps 1-5 for the next Cal.
- Mix Calibrator well by gently shaking for 1 minute.

Refer to the Stat Profile Prime Plus VET Analyzer Instructions for Use Manual for complete directions. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents. Cartridge is ready for use.

Traceability of Standards

Analytes traced to NIST Standard Reference Materials wherever possible.

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