REF 52427, 52861, 53365, 53105, 53469, 53470, 52862, 52863, 53364, 53466, 53467, 53468

Rx Only IVD 15°C - 30°C CE

# Stat Profile® Prime Calibrator Cartridge CCS/CCS Comp

Kalibratorkassette, Κασέτα Βαθμονομητή, Cartucho del calibrador, Cartouche étalon, Cartuccia calibratore, Calibrator Cartridge, Kalibratorkassett, Kalibrátorkazetta, キャリブレータ カートリッジ, 측정기 카트리지, 校准液盒

LOT 23198016

2024-10-27

CAL A	≥ 500 mL	CAL B	≥ 250 mL	CAL F	≥ 720 mL	SOLN R	≥ 300 mL
pH	7.349	pH	6.840	PCO <sub>2</sub>	51.3 mmHg	KCl	2 mol/L
PCO <sub>2</sub>	25.1 mmHg	PO <sub>2</sub>	167.0 mmHg	PO <sub>2</sub>	93.5 mmHg	 WASTE CONTAINER W	
Na <sup>+</sup>	133.0 mmol/L	Na <sup>+</sup>	68.5 mmol/L				
K <sup>+</sup>	4.0 mmol/L	K <sup>+</sup>	10.0 mmol/L				
Cl <sup>-</sup>	106.0 mmol/L	Cl <sup>-</sup>	46.0 mmol/L				
iCa	1.10 mmol/L	iCa	2.20 mmol/L				
Glu	80 mg/dL	Glu	200 mg/dL				
Lac	2.0 mmol/L	Lac	10.0 mmol/L				

**EN**  
**Product Description**  
 The above table lists the concentrations of constituents for the calibrators and solutions contained in the Calibrator Cartridge. Contains the following aqueous solutions within individual flexible bags housed in a cardboard box: 3 Calibrators (Calibrator A, B, and F), reference solution, and a flexible waste bag. Each bag includes a firmet with a septa. The exposed bag firmets are arranged in a line along the rear of the cartridge box. These septa are pierced during the insertion of the cartridge into the analyzer.

**Intended Use**  
 The Stat Profile Prime Calibrator Cartridge CCS is intended for the calibration of pH, PCO<sub>2</sub>, PO<sub>2</sub>, Hct, Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, iCa, Glucose (Glu), and Lactate (Lac) using the Stat Profile Prime CCS Analyzer.

**Methodology**  
 Refer to Prime Analyzer Instructions For Use Manual for Test Methodologies and Principles.

**Composition**  
 Calibrator volumes are listed in the table above. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF: NCCLS DOCUMENT M29-T2).

**Warnings and Cautions:**  
 Intended for in vitro diagnostic use.  
 Refer to the Prime Analyzer Instructions for Use Manual for complete directions for use, installation instructions, recommendations for calibration intervals, troubleshooting information, Methodology and Principles of the testing procedures. Follow standard practices required for handling laboratory reagents. Discard Calibrator Cartridges and Waste products in accordance with local guide lines. Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

**Storage**  
 Store at 15-30°C. Expiration date is printed on the front of the Calibrator Cartridge. Directions For Use Refer to the Prime Analyzer Instructions For Use for complete directions for installation of the cartridge, and recommendations for Calibration Intervals. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents. Cartridge is ready for use.

**Traceability of Standards**  
 Analytes traced to NIST Standard Reference Materials.

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