REF 57825, 57826, 57827, 57828, 57829, 65141

Rx only IVD 15°C 30°C CE

Stat Profile Prime Plus® Calibrator Cartridge

Stat Profile Prime Plus® Kalibrierkassette, Φυσιγγιό μέσου βαθμονόμησης με κρεατινίνη Stat Profile Prime Plus®, Cartucho de calibrador Stat Profile Prime Plus®, Cartouche de calibrateurs Stat Profile Prime Plus®, Cartuccia calibratore Stat Profile Prime Plus®, Cartucho de calibrador Stat Profile Prime Plus®, Stat Profile Prime Plus® kalibrator patron, ייְהִוָּתְהָרָה Stat Profile Prime Plus®, Stat Profile Prime Plus® カリブレーターカートリッジ, Stat Profile Prime Plus® 보정기 카트리지, Stat Profile Prime Plus® 校准试剂盒

LOT 23209042

2025-01-13

CAL A	CAL B	CAL C	CAL D	CAL F	SOLN R				
HCO ₃ Na ⁺ K ⁺ Cl ⁻ iCa iCa iMg iMg Glu Glu Lac Lac BUN BUN	11 mmol/L 135.5 mmol/L 4.0 mmol/L 108.0 mmol/L 1.04 mmol/L 4.16 mg/dL 0.50 mmol/L 1.22 mg/dL 80 mg/dL 4.44 mmol/L 2.0 mmol/L 17.8 mg/dL 10 mg/dL 3.6 mmol/L	pH H ⁺ Na ⁺ K ⁺ iCa iCa iMg iMg tHb tHb tHb tHb	6.834 @ 37°C 146.55 nmol/L 77.0 mmol/L 10.0 mmol/L 1.92 mmol/L 7.68 mg/dL 1.50 mmol/L 3.63 mg/dL 15.0 g/dL 150 g/L 9.3 mmol/L	HCO ₃ Cl ⁻ Glu Glu Lac Lac BUN BUN	22 mmol/L 69.0 mmol/L 200 mg/dL 11.10 mmol/L 10.0 mmol/L 89.0 mg/dL 50 mg/dL 17.9 mmol/L	pH H ⁺ PO ₂ PO ₂	7.354 @ 37°C 44.26 nmol/L 102.5 mmHg 13.63 kPa	KCl WASTE CONTAINER W	2 mol/L

EN

Description

The Cartridge contains the following aqueous solutions within individual bags: 5 Calibrators (Calibrator A, B, C, D, and F SOLN), Reference Solution, and a flexible waste bag. Each bag includes a filament with a septum. The exposed bag filaments are arranged in a line along the rear of the cartridge box. The septa are pierced during the insertion of the cartridge into the analyzer.

Intended Use

For the quantitative determination of Na⁺, K⁺, Cl⁻, iCa, iMg, BUN (Urea), Glucose, Lactate, pH, PCO₂, PO₂, SO₄²⁻, Hct, Hb, Hbf*, O₂Hb, COHb, MetHb, HHb, and iBil* in human blood using the Stat Profile Prime Plus Analyzer.

*Not available in the USA or for Point-of-Care/Near-Patient Testing use.

Methodology

Refer to Stat Profile Prime Plus Analyzer Instructions For Use Manual for Methodologies and Principles.

Composition

The cartridges are aqueous salt solutions containing analytes at concentrations listed in the table above. Calibrator B also contains dye. Contain no constituents of human origin, however, good laboratory practice should be followed during the handling of these materials. (REF. NCLS DOCUMENT M29-T2).

Intended for In-vitro diagnostic use

Refer to the Stat Profile Prime Plus Analyzer Instructions For Use Manual for complete directions. Follow standard practices required for handling laboratory reagents. Discard Cartridge in accordance with local guidelines.

Vial Bag Contains Preservative (eye, skin, and respiratory irritant).

Storage

Store at 15 - 30°C. Expiration date is printed on the Calibrator Cartridge.

Directions For Use

Refer to the Stat Profile Prime Plus Analyzer Instructions For Use for complete directions. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents.

Traceability of Standards

Analytes traced to NIST Standard Reference Materials wherever possible.

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