REF 57831, 57832, 57833, 57834, 57835, 65142

Rx only IVD 15°C -30°C CE

Stat Profile Prime Plus® Calibrator Cartridge With Creatinine

Stat Profile Prime Plus® Kalibrierkassette mit Kreatinin, Stat Profile Prime Plus® Φυσιγγίου μέσου βαθμονόμησης με κρεατινίνη, Cartucho de calibrador con creatinina Stat Profile Prime Plus®, Cartouche de calibrateurs Stat Profile Prime Plus® avec créatinine, Cartuccia calibratore con creatinina Stat Profile Prime Plus®, Cartucho de calibrador Stat Profile Prime Plus® com creatinina, Stat Profile Prime Plus® kalibrátor patron kreatininnel, ژیستاپرایم پرم جم مهندسی Stat Profile Prime Plus®, Stat Profile Prime Plus® キャリブレーターカートリッジ (creatine), 콜레아티닌 사용 Stat Profile Prime Plus® 보정기 카트리지, Stat Profile Prime Plus® 含肌酐的校准试剂盒

LOT 23209041

2025-01-11

CAL A	CAL B	CAL C	CAL D	CAL F	SOLN R
HCO ₃	11 mmol/L	pH 6.834 @ 37°C	HCO ₃	22 mmol/L	pH 7.354 @ 37°C
Na ⁺	135.5 mmol/L	H ⁺ 146.55 nmol/L	Cl ⁻	69.0 mmol/L	H ⁺ 44.26 nmol/L
K ⁺	4.0 mmol/L	Na ⁺ 77.0 mmol/L	Glu	200 mg/dL	PO ₂ 102.5 mmHg
Cl ⁻	108.0 mmol/L	K ⁺ 10.0 mmol/L	Glu	11.10 mmol/L	PO ₂ 13.63 kPa
iCa	1.04 mmol/L	iCa	Lac	10.0 mmol/L	KCl 2 mol/L
iCa	4.16 mg/dL	iCa	Lac	89.0 mg/dL	
iMg	0.50 mmol/L	iMg	BUN	50 mg/dL	
iMg	1.22 mg/dL	iMg	BUN	17.9 mmol/L	
Glu	80 mg/dL	tHb	Creat	5.0 mg/dL	
Glu	4.44 mmol/L	tHb	Creat	0.44 mmol/L	
Lac	2.0 mmol/L	tHb	Creat	442 μmol/L	
Lac	17.8 mg/dL				
BUN	10 mg/dL				
BUN	3.6 mmol/L				
Creat	1.0 mg/dL				
Creat	0.09 mmol/L				
Creat	88.4 μmol/L				

EN**Description**

The cartridge contains the following aqueous solutions within individual bags: 5 Calibrators (Calibrator A, B, C, D, and F), Reference Solution, and a flexible waste bag. Each bag includes a filament with a septum. The exposed bag filaments are arranged in a line along the rear of the cartridge box. The exposed bag filaments are pierced during the insertion of the cartridge into the analyzer. A charge kit consisting of 2 syringes filled with creatione solution is included for cartridge activation prior to installation.

Intended Use

For the quantitative determination of Na⁺, K⁺, Cl⁻, iCa, iMg, BUN (Urea), Glucose, Lactate, Creatinine, pH, PCO₂, PO₂, SO₄²⁻, Hct, Hb, HbF, O₂Hb, COHb, MetHb, HHb, and tBil¹ in human blood using the Stat Profile Prime Plus Analyzer.

Methodology

Refer to Stat Profile Prime Plus Analyzer Instructions for Use Manual for Methodologies and Principles.

Compatibility

The calibrators are aqueous salt solutions containing analytes at concentrations listed in the table above. Calibrator B also contains dyes. Contain no constituents of human origin; however, good laboratory practice should be followed during the handling of these materials. (REF. NCCLS DOCUMENT M29-T2).

Warnings and Caution

Intended for *in vitro* diagnostic use.
Refer to the Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory reagents. Discard Cartridge in accordance with local guidelines.

Waste Bag Contains Preservative (e.g., skin, and respiratory irritant).

Storage: Store at 15–30°C. Expiration date is printed on the Calibrator Cartridge.

Directions For Use**Action Items for Use**

Cal A and Cal C must be charged with creatinine syringes prior to installation of the Calibrator Cartridge. The two syringes are labeled and color coded to correspond to the filaments for the cartridge for Cal A and Cal C.

1. Turn the syringe with tip side down, and remove protective cap.
2. Attach one of the standard needle assemblies to the syringe. Remove the protective cover from the needle.
3. Match the color and label of the syringe to the appropriate filament and insert needle.
4. Slightly depress syringe plunger until the contents are dispensed.
5. DO NOT PULL BACK ON THE PLUNGER TO FLUSH CONTENTS OF SYRINGE.
6. Remove needle/syringe assembly from filament and discard in an appropriate sharps container.
7. Repeat Steps 1-5 for the next Cal.
7. Mix Cartridge well by gently shaking for 1 minute.

Refer to the Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete directions. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents. Cartridge is ready for use.

*Not available in the USA or for Point-of-Care/Near-Patient Testing use.

Traceability of Standards

Analytes traced to NIST Standard Reference Materials wherever possible.

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