REF 53104, 53360, 53359, 53463, 53464, 53465

Rx IVD 15°C 30°C CE

Stat Profile® Prime Calibrator Cartridge ABC

Kalibratorkassette, Κασέτα Βαθμονομητή, Cartucho del calibrador, Cartucho étonal, Cartuccia calibratore, Calibrator Cartridge, Kalibratorkassett

LOT 23229031

2025-02-02

CAL A ≥ 500 mL	CAL B ≥ 250 mL	CAL F ≥ 720 mL	SOLN R ≥ 300 mL
pH 7.350	pH 6.840	PCO ₂ 51.6 mmHg	KCl 2 mol/L
PCO ₂ 25.2 mmHg	PO ₂ 167.0 mmHg	PO ₂ 93.8 mmHg	WASTE CONTAINER W

Product Description
The above table lists the concentrations of constituents for the calibrators and solutions contained in the Calibrator Cartridge. Contains the following aqueous solutions within individual flexible bags housed in a cartridge box: 3 Calibrators (Calibrator A, B, and F), Reference solution, and a flexible waste bag. Each bag includes a fitting with a septa. The expandable bag is easily arranged in a line along the rest of the cartridge box. These septa are pierced during the insertion of the cartridges into the analyzer.

Intended Use
For the quantitative determination of pH, PCO₂, and PO₂ in human blood using the Nova Biomedical Prime ABC Analyzer.

Methodology
Refer to Prime Analyzer Instructions For Use Manual for Test Methodologies and Principles.

Composition
Calibrator volumes are listed in the table above. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF. NCCLS DOCUMENT M29-T2).

Warnings and Cautions:
Intended for In vitro diagnostic use.
Refer to Prime Analyzer Instructions for Use Manual for complete directions for use, installation instructions, recommendations for calibration intervals, troubleshooting information, Methodology and Principles of the testing procedures. Follow standard practices required for handling laboratory reagents. Discard Calibrator Cartridges and Waste products in accordance with local guidelines. Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

Storage
Store at 15-30°C. Expiration date is printed on the front of the Calibrator Cartridge.
Directions For Use Refer to the Prime Analyzer Instructions for Use for complete directions for installation of the cartridge, and recommendations for Calibration Intervals. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents. Cartridge is ready for use.

Traceability of Standards
Analytes traced to NIST Standard Reference Materials.

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