REF 58404, 58405

15°C | 30°C CE

Stat Profile Prime Plus® VET Calibrator Cartridge With Creatinine

Stat Profile Prime Plus® VET Kalibrierkassette mit Kreatinin, Φυσιγγιο μέσου βαθμονόμησης VET με κρεατίνην Stat Profile Prime Plus®, Cartucho de calibrador con creatinina Stat Profile Prime Plus® VET, Cartouche de calibrateurs avec créatinine Stat Profile Prime Plus® VET, Cartuccia calibratore con creatinina Stat Profile Prime Plus® VET, Cartucho de calibrador Stat Profile Prime Plus® VET com creatinina, Stat Profile Prime Plus® VET kalibrátor patron kreatininom, Stat Profile Prime Plus® VET קרטינין כוונתית ביל' עט קרייטינין צייל בערבית, Stat Profile Prime Plus® VET キャリブレータ(クリアチニン), 크레아티닌 사용 Stat Profile Prime Plus® VET 보정기 카트리지, Stat Profile Prime Plus® 含肌酐的校准试剂盒

LOT 23263040

2025-02-07

| CAL A | ≥ 495 mL | CAL B | ≥ 250 mL | CAL C | ≥ 245 mL | CAL D | ≥ 250 mL | CAL F | ≥ 850 mL | SOLN R | ≥ 500 mL |
|-------------------------------|--------------|-----------------|---------------|-------------------------------|--------------|----------------|--------------|-----------------|------------|----------|-------------------|
| HCO ₃ ⁻ | 11 mmol/L | pH | 6.834 @ 37°C | HCO ₃ ⁻ | 22 mmol/L | pH | 7.354 @ 37°C | PO ₂ | 104.7 mmHg | KCl | 2 mol/L |
| Na ⁺ | 135.5 mmol/L | H ⁺ | 146.55 nmol/L | Cl ⁻ | 69.0 mmol/L | H ⁺ | 44.26 nmol/L | PO ₂ | 13.93 kPa | | |
| K ⁺ | 4.0 mmol/L | Na ⁺ | 77.0 mmol/L | Glu | 200 mg/dL | | | | | ! | WASTE CONTAINER W |
| Cl ⁻ | 108.0 mmol/L | K ⁺ | 10.0 mmol/L | Glu | 11.10 mmol/L | | | | | | |
| iCa | 1.04 mmol/L | iCa | 1.92 mmol/L | Lac | 10.0 mmol/L | | | | | | |
| iCa | 4.16 mg/dL | iCa | 7.68 mg/dL | Lac | 89.0 mg/dL | | | | | | |
| IMg | 0.50 mmol/L | IMg | 1.50 mmol/L | BUN | 50 mg/dL | | | | | | |
| IMg | 1.22 mg/dL | IMg | 3.63 mg/dL | BUN | 17.9 mmol/L | | | | | | |
| Glu | 80 mg/dL | tHb | 15.4 g/dL | Creat | 5.0 mg/dL | | | | | | |
| Glu | 4.44 mmol/L | tHb | 154 g/L | Creat | 0.44 mmol/L | | | | | | |
| Lac | 2.0 mmol/L | tHb | 9.6 mmol/L | Creat | 442 µmol/L | | | | | | |
| Lac | 17.8 mg/dL | | | | | | | | | | |
| BUN | 10 mg/dL | | | | | | | | | | |
| BUN | 3.6 mmol/L | | | | | | | | | | |
| Creat | 1.0 mg/dL | | | | | | | | | | |
| Creat | 0.09 mmol/L | | | | | | | | | | |
| Creat | 88.4 µmol/L | | | | | | | | | | |

Description

The Cartridge contains the following aqueous solutions within individual bags: 5 Calibrators (Calibrator A, B, C, D, and F), Reference Solution, and a flexible waste bag. Each bag includes a unique color and code to correspond to the fittings on the cartridge and the rest of the cartridge box. The septa are pierced during the insertion of the cartridge into the analyzer. A charge kit consisting of 2 syringes filled with creatinine solution is included for calibrator activation prior to installation.

Intended Use

For the quantitative determination of Na⁺, K⁺, Cl⁻, iCa, IMg, BUN (Urea), Glucose, Lactate, Creatinine, pH, PO₂, PO₃, SO₄, Hct, tHb, COHb, MetHb, HHb, and tBil in animal whole blood using the Stat Profile Prime Plus VET Analyzer.

Methodology

Refer to Stat Profile Prime Plus VET Analyzer Instructions For Use Manual for Methodologies and Principles.

Composition

Calibrator and solution volumes are listed in the table. Contains no constituents of human origin; however, good laboratory practice should be followed during handling of these materials. (REF NCCLS DOCUMENT M26-T2).

Warnings and Cautions

Refer to the Stat Profile Prime Plus VET Analyzer Instructions For Use Manual for complete directions. Do not mix with other reagents or handling laboratory reagents. Discard Cartridge in accordance with local guidelines.

Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

Storage

Store at 15 - 30°C. Expiration date is printed on the Calibrator Cartridge.

Directions For Use

Activating the Cartridge

Cal A and Cal C must be charged with creatinine syringe prior to installation of the Calibrator Cartridge. The two syringes are labeled and color coded to correspond to the fittings on the carton for Cal A and Cal C.

1. Remove syringe with tip side down and remove protective cap.
2. Attach the one enclosed needle assembly to the syringe. Remove the protective cover from the needle.
3. Match the color and label of the syringe to the appropriate fitting and insert needle.
4. Slowly depress syringe plunger until the contents are dispensed.
5. DO NOT PULL BACK ON THE PLUNGER TO FLUSH CONTENTS OF SYRINGE.
6. Repeat needle/syringe assembly from fitting and discard in an appropriate sharps container.
7. Mix Cartridge well by gently shaking for 1 minute.

Refer to the Stat Profile Prime Plus VET Analyzer Instructions For Use Manual for complete directions. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents. Cartridge is ready for use.

Traceability of Standards

Analyses traced to NIST Standard Reference Materials wherever possible.

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Waltham, MA 02454 U.S.A.