REF 58395,58396

15°C 30°C CE

# Stat Profile Prime Plus® VET Calibrator Cartridge

Stat Profile Prime Plus® VET Kalibrierkassette, Φυσιγγίο βαθμονόμησης VET Stat Profile Prime Plus®, Cartucho de calibrador Stat Profile Prime Plus® VET, Cartouche de calibrateur Stat Profile Prime Plus®, Stat Profile Prime Plus® VET kalibrátor patron, Stat Profile Prime Plus® VET, Stat Profile Prime Plus® VET, Stat Profile Prime Plus® VET キャリブレータートリッジ, Stat Profile Prime Plus® VET 보정기 카트리지, Stat Profile Prime Plus® VET 校准试剂盒

LOT 23270044

2025-03-14

CAL A	$\geq 495 \text{ mL}$	CAL B	$\geq 250 \text{ mL}$	CAL C	$\geq 245 \text{ mL}$	CAL D	$\geq 250 \text{ mL}$	CAL F	$\geq 850 \text{ mL}$	SOLN R	$\geq 500 \text{ mL}$
HCO <sub>3</sub>	11 mmol/L	pH	6.834 @ 37°C	HCO <sub>3</sub>	22 mmol/L	pH	7.354 @ 37°C	PO <sub>2</sub>	106.6 mmHg	KCl	2 mol/L
Na <sup>+</sup>	135.5 mmol/L	H <sup>+</sup>	146.55 nmol/L	Cl <sup>-</sup>	69.0 mmol/L	H <sup>+</sup>	44.26 nmol/L	PO <sub>2</sub>	14.18 kPa		WASTE CONTAINER W
K <sup>+</sup>	4.0 mmol/L	Na <sup>+</sup>	77.0 mmol/L	Glu	200 mg/dL						
Cl <sup>-</sup>	108.0 mmol/L	K <sup>+</sup>	10.0 mmol/L	Glu	11.10 mmol/L						
iCa	1.04 mmol/L	iCa	1.92 mmol/L	Lac	10.0 mmol/L						
iCa	4.16 mg/dL	iCa	7.68 mg/dL	Lac	89.0 mg/dL						
iMg	0.50 mmol/L	iMg	1.50 mmol/L	BUN	50 mg/dL						
iMg	1.22 mg/dL	iMg	3.63 mg/dL	BUN	17.9 mmol/L						
Glu	80 mg/dL	tHb	15.4 g/dL								
Glu	4.44 mmol/L	tHb	154 g/L								
Lac	2.0 mmol/L	tHb	9.6 mmol/L								
Lac	17.8 mg/dL										
BUN	10 mg/dL										
BUN	3.6 mmol/L										

## EN

## Description

The Cartridge contains the following aqueous solutions within individual bags: 5 Calibrators (Calibrators A, B, C, D, and F), Reference Solution, and a flexible waste bag. Each bag includes a unique identification number. The enclosed bag fittings are arranged in a line along the rear of the cartridge box. These sets are placed during the insertion of the cartridge into the analyzer. A charge kit consisting of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation.

## Intended Use

For the quantitative determination of Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, iCa, Mg, BUN (Urea), Glucose, Lactate, pH, PO<sub>2</sub>, PO<sub>3</sub>, SO<sub>4</sub>%, Hct, tHb, O<sub>2</sub>Hb, COHb, MetHb, HHb, and iMg in animal whole blood using the Stat Profile Prime Plus VET Analyzer.

## Methodology

Refer to Stat Profile Prime Plus VET Analyzer Instructions For Use Manual for Methodologies and Principles.

## Composition

Calibrator and solution volumes are listed in the table. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF: NCCLS DOCUMENT M29-T2).

## Warnings and Cautions

Refer to the Stat Profile Prime Plus VET Analyzer Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory reagents. Discard cartridge in accordance with local guidelines.

Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

## Storage

Store at 15 - 30°C. Expiration date is printed on the Calibrator Cartridge.

## Directions For Use

Refer to the Stat Profile Prime Plus VET Analyzer Instructions for Use for complete directions. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents.

## Traceability of Standards

Analytes traced to NIST Standard Reference Materials wherever possible.

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