REF 57831, 57832, 57833, 57834, 57835, 65142

Rx only IVD 15°C -30°C CE

# Stat Profile Prime Plus® Calibrator Cartridge With Creatinine

Stat Profile Prime Plus® Kalibrierkassette mit Kreatinin, Stat Profile Prime Plus® Φυσιγγίου μέσου βαθμονόμησης με κρεατινίνη, Cartucho de calibrador con creatinina Stat Profile Prime Plus®, Cartouche de calibrateurs Stat Profile Prime Plus® avec créatinine, Cartuccia calibratore con creatinina Stat Profile Prime Plus®, Cartucho de calibrador Stat Profile Prime Plus® com creatinina, Stat Profile Prime Plus® kalibrátor patron kreatininnel, קרייניןן מודול Stat Profile Prime Plus®, Stat Profile Prime Plus® キャリブレーターカートリッジ (ケレアチニン), 크레아티닌 사용 Stat Profile Prime Plus® 보정기 카트리지, Stat Profile Prime Plus® 含肌酐的校准试剂盒

LOT 23283017

2025-03-22

CAL A	CAL B	CAL C	CAL D	CAL F	SOLN R
HCO <sub>3</sub>	11 mmol/L	pH 6.834 @ 37°C	HCO <sub>3</sub> 22 mmol/L	pH 7.354 @ 37°C	PO <sub>2</sub> 104.6 mmHg
Na <sup>+</sup>	135.5 mmol/L	H <sup>+</sup> 146.55 nmol/L	Cl <sup>-</sup> 69.0 mmol/L	H <sup>+</sup> 44.26 nmol/L	PO <sub>2</sub> 13.91 kPa
K <sup>+</sup>	4.0 mmol/L	Na <sup>+</sup> 77.0 mmol/L	Glu 200 mg/dL		
Cl <sup>-</sup>	108.0 mmol/L	K <sup>+</sup> 10.0 mmol/L	Glu 11.10 mmol/L		
iCa	1.04 mmol/L	iCa 1.92 mmol/L	Lac 10.0 mmol/L		
iCa	4.16 mg/dL	iCa 7.68 mg/dL	Lac 89.0 mg/dL		
iMg	0.50 mmol/L	iMg 1.50 mmol/L	BUN 50 mg/dL		
iMg	1.22 mg/dL	iMg 3.63 mg/dL	BUN 17.9 mmol/L		
Glu	80 mg/dL	tHb 15.5 g/dL	Creat 5.0 mg/dL		
Glu	4.44 mmol/L	tHb 155 g/L	Creat 0.44 mmol/L		
Lac	2.0 mmol/L	tHb 9.6 mmol/L	Creat 442 µmol/L		
Lac	17.8 mg/dL				
BUN	10 mg/dL				
BUN	3.6 mmol/L				
Creat	1.0 mg/dL				
Creat	0.09 mmol/L				
Creat	88.4 µmol/L				

**EN**

**Description**  
The Cartridge contains the following aqueous solutions within individual bags: 5 Calibrators (Calibrator A, B, C, D, and F), Reference water, and 10ml preservative. The bag includes a septum and a needle. The exposed bag fibres are arranged in a line along the rear of the cartridge body. The septa are pierced during the insertion of the cartridge into the analyzer. A charge kit consisting of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation.

**Intended Use**

For the quantitative determination of Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, iCa, iMg, BUN (Urea), Glucose, Lactate, Creatinine, pH, PCO<sub>2</sub>, PO<sub>2</sub>, SO<sub>4</sub><sup>2-</sup>, Hct, tHb, HbF, O<sub>2</sub>Hb, COHb, MetHb, HHb, and iBil<sup>+</sup> in human blood using the Stat Profile Prime Plus Analyzer.

**Methodology**

Refer to Stat Profile Prime Plus Analyzer Instructions for Use Manual for Methodologies and Principles.

**Composition**

The calibrators are aqueous salt solutions containing analytes at concentrations listed in the table above. Calibrator B also contains dye. Contain no constituents of human origin, however laboratory practice should be followed during the handling of these materials. (REF. NCCLS DOCUMENT M29-2).

**Warnings and Precautions**

Indicated for in vitro diagnostic use.  
Refer to the Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete directions. Standard practices required for handling laboratory reagents. Discard Cartridge in accordance with local guidelines.

Vaste Bag Contains: Preservative (eye, skin, and respiratory irritant).

**Storage**

Store at 15 - 30 °C. Expiration date is printed on the Calibrator Cartridge.

**Directions For Use****Activation of Cartridge**

Cal A and Cal C must be charged with creatinine syringe prior to installation of the Calibrator Cartridge. The two syringes are labeled and color coded to correspond to the fittings on the cartridge body.

1. Hold the syringe with tip side down and remove protective cap.

2. Attach one of the enclosed needle assemblies to the syringe. Remove the protective cover from the needle.

3. Slowly depress the label of the syringe to the appropriate fitting and insert needle.

4. Slowly depress syringe plunger until the contents are dispensed.

DO NOT PULL BACK ON THE PLUNGER TO FLUSH CONTENTS OF SYRINGE.

5. Remove needle/syringe assembly from fitting and discard in an appropriate sharps disposal container.

6. Repeat Steps 1-5 for the next Cal.

7. Mix cartridge well by gently shaking for 1 minute.

Refer to the Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete directions. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents. Cartridge is ready for use.

\*Not available in the USA or for Point-of-Care/Near-Patient Testing use.

**Traceability of Standards**

Analyses traced to NIST Standard Reference Materials wherever possible.

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**EL**

**Properties**  
The Cartridge contains the following aqueous solutions in individual bags: 5 Calibrators (Calibrator A, B, C, D and F), Reference water, and 10ml desinfectant. The bag includes a septum and a needle. The exposed bag fibers are arranged in a line along the rear of the cartridge body. The septa are pierced during the insertion of the cartridge into the analyzer. A charge kit consisting of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation.

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