Stat Profile® pHox Ultra/Critical Care Xpress Blood Gas Controls Auto-Cartridge

Blutgas-Kontroll-Auto-Kassetten, Αυτόματη Κασέτα Διαλυμάτων Ελέγχου Αερίων Αίματος, Cartucho automático para controles de gases en sangre, Cartouche CQ pour les gaz du sang, Cartuccia automatica controlli, Cartucho automático de controle de gases no sangue, Blodgaskontroller, autokassett

CONTROL 1 2 3

LOT 23271019 2025-03-11

Table with 4 columns: Parameter, Unit, CONTROL 1 (min - x̄ - max), CONTROL 2 (min - x̄ - max), CONTROL 3 (min - x̄ - max). Rows include pH, H+, pCO2, pCO2, pO2, pO2, SO2, Hct, Hb, Hb.

Product Description An aqueous quality control material for monitoring the measurement of pH, pCO2, pO2, SO2, hematoctrit (Hct), and hemoglobin (Hb) for use with Nova Biomedical analyzers ONLY.

Acidosis, with low SO2, low-normal Hct/Hb Normal pH Alkalosis, with SO2, High Hct/Hb

Intended Use For in vitro diagnostic use for monitoring the performance of Nova Biomedical Stat Profile pHox Ultra/Critical Care Xpress Analyzers.

Methodology Refer to Stat Profile pHox Ultra/CX Analyzer Instructions For Use Manual.

Composition A buffered bicarbonate solution, the Controls Auto-Cartridge consists of 3 flexible bags, each with a known pH. Solutions are equilibrated with known levels of O2, CO2, and N2.

Warnings and Cautions: Must be stored at 2-8°C for at least 24 hours before using. DO NOT SHAKE the cartridge.

Storage Store at 2-8°C; DO NOT FREEZE. It is critical to follow the room temperature equilibrium instructions prior to use as described in "Directions for Use".

Directions for use Controls must be stored at approximately 24-26°C for at least 24 hours prior to opening.

Limitations PO2 values vary inversely with temperature (approximately 1%°C).

Traceability of Standards Analytes are traced to NIST Standard Reference Materials.

Reference Intervals Concentrations are formulated to represent three pH levels (Acidosis, Normal pH, and Alkalosis), and with Low-Normal and High Hct/Hb levels.

Expected Ranges The EXPECTED RANGE for each analyte was determined at Nova Biomedical by using multiple runs of each level of control at 37°C on multiple instruments.

How to Define and Determine Reference Intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13

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