Stat Profile® Prime Auto QC Cartridge ABG

Kontroll-Auto-Kassetten, Autómatu Kaoščta Διαλυμάτων Ελέχου, Cartucho automático para controles, Cartouche CQ, Cartuccia automatica controlli, Cartucho automático de controlo, Kontroller, autokassett

CONTROL | 1 | 2 | 3

LOT 23262019 2025-02-21

Expected Ranges, Erwartete Bereiche, Αναμενόμενα Εύρη, Rangos esperados, Fourchettes attendues, Intervalli previsti, Gamas previstas, Förväntade områden				
		CONTROL 1 min - \bar{x} - max	CONTROL 2 min - \bar{x} - max	CONTROL 3 min - \bar{x} - max
pH		7.138 - 7.168 - 7.198	7.325 - 7.355 - 7.385	7.560 - 7.590 - 7.620
H ⁺	nmol/L	72.78 - 67.92 - 63.39	47.32 - 44.16 - 41.21	27.54 - 25.70 - 23.99
pCO ₂	mmHg	58.0 - 65.0 - 72.0	41.0 - 46.0 - 51.0	18.8 - 22.8 - 26.8
pCO ₂	kPa	7.73 - 8.67 - 9.60	5.47 - 6.13 - 6.80	2.51 - 3.04 - 3.57
pO ₂	mmHg	46.4 - 56.4 - 66.4	85.9 - 95.9 - 105.9	126.5 - 141.5 - 156.5
pO ₂	kPa	6.19 - 7.52 - 8.85	11.45 - 12.79 - 14.12	16.87 - 18.87 - 20.86

Product Description
Consists of 3 flexible bags within a cardboard carton. Each bag contains an aqueous quality control material for monitoring the measurement of pH, PCO₂, and PO₂ for use with Nova Biomedical analyzers ONLY.

Formulated at three levels:

CONTROL | Acidosis

CONTROL | Normal pH

CONTROL | Alkalosis

Intended Use
The Stat Profile Prime Auto QC Cartridge ABG is a Quality Control Material intended for *in vitro* diagnostic testing by healthcare professionals for monitoring the performance of the Stat Profile Prime ABG.

Methodology
Refer to Prime Analyzer Instructions For Use Manual.

Composition

A buffered bicarbonate solution, each control with a known pH. Solutions are equilibrated with known levels of CO₂ and O₂. Measured. Each bag contains a minimum volume of 100 mL. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF. NCCLS DOCUMENT M29-T2).

Warnings and Cautions:

Material must be stored at 24°C for at least 24 hours before use. DO NOT FREEZE. Intended for *in vitro* diagnostic use. Refer to the Prime Analyzer Instructions for Use for complete directions for use, including recommendations for use of controls, troubleshooting information, Methodology and Principles of the testing procedure. Follow standard practices required for handling laboratory reagents.

Storage:

Store at 0-8°C; DO NOT FREEZE. Expiration date is printed on each cartridge.

Directions for use

Cartridge must be stored at approximately 24-26°C for at least 24 hours prior to opening. Refer to Analyzer Instructions for Use for complete instructions. Verify that the lot number appearing on the Expected Ranges Table is this same lot number indicated on the control cartridge. The cartridge should be stored by gently inverting for several seconds. Do not shake the cartridge.

Limitations:

PO₂ values vary inversely with temperature (approximately 1%/°C).

The Expected Range values are specific for instruments and calibrators manufactured by Nova Biomedicals.

Traceability of Standards

Analytes are traced to NIST Standard Reference Materials.

Reference Intervals

Concentrations are formulated to represent true pH levels (Acidosis, Normal pH, and Alkalosis). The expected clinical range value for these analytes in patient blood is referenced in Tietz, NW ed 1988 Textbook of Clinical Chemistry, W.B. Saunders Co.

Users may wish to determine MEAN VALUES and EXPECTED RANGES in their own laboratory.

Expected Range:

The EXPECTED RANGE for each analysis was determined at Nova Biomedical

by use of the means of each level of control at 37°C on multiple instruments.

The EXPECTED RANGE indicates the maximum deviations from the mean value that may be expected under differing laboratory conditions for instruments operating within specifications.

Refer to Expected Ranges Table.

"How to Define and Determine Reference Intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13

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