Stat Profile® pHox Ultra/Critical Care Xpress Chemistry Controls & Performance Check Solution Auto-Cartridge With Creatinine

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Product Description
Carton contains three bags. Two bags contain assayed quality control material formulated for monitoring the measurement of Na⁺, K⁺, Cl⁻, Ca⁺⁺, Mg⁺⁺, Glucose, Lactate, BUN and Creatinine. One bag contains a buffered salt solution for glucose and creatinine membrane performance (Performance Check Solution). A change kit consisting of three syringes filled with creatinine solution is included for cartridge activation prior to installation.

Intended use
For in vitro diagnostic use to monitor the performance of the Stat Profile pHox Ultra/Critical Care Xpress Analyzer and to insure optimum performance of glucose and creatinine membranes. Performance check solution should be analyzed before controls are run and after changing glucose or creatinine membranes.

Methodology
Refer to Stat Profile pHox Ultra/Critical Care Xpress Analyzer Instructions for Use Manual.

Composition
Controls (4 and 5) contain known concentrations of Na⁺, K⁺, Cl⁻, Ca⁺⁺, Mg⁺⁺, Glucose, Lactate, BUN, and Creatinine at two levels. Performance Check Solution is a buffered salt solution containing reagents for glucose and creatinine membrane analysis. Each bag contains a minimum volume of 100 mL.

Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials.

Warning and Precautions
For use on Nova Biomedical Analyzers. Do not freeze. Refer to the Stat Profile pHox Ultra/Critical Care Xpress Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory controls.

Storage
Store at 2-8°C. Do not freeze. Each cartridge has lot number and expiration date printed on the label.

Directions for Use

Action of Controls
Each level of control must be charged with the creatinine syringe prior to installation of the Controls Auto-Cartridge. The three syringes are labeled and color coded to correspond to the three levels of control for each cartridge.

1. Hold the syringe with tip side down and remove protective cap.
2. Attach one of the enclosed needle assemblies into the syringe. Remove the protective cover from the needle assembly.
3. Match the color and label of the syringe to the appropriate fitment and insert needle.

4. Slowly pull back on the plunger until the desired dispensing.
5. Remove needle/syringe assembly from fitment and discard in appropriate sharps container.
6. Repeat Steps 1-5 for next two levels.
7. Mix the contents of the cartridge for 1 minute.
DO NOT SHAKE THE CARTRIDGE.

Limitations
For use on Nova Biomedical Analyzers Only. After installation, the Auto-Cartridge may be used for a maximum of 21 days from the date of initial activation and installation. The Analyzer will indicate that the cartridge is invalid at that time.

This product is warranted to perform as described in this product insert, and Nova Biomedical disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall Nova Biomedical be liable for any consequential damages arising out of the aforesaid expressed warranty.

Traceability of Controls/Solution

Analyses are traceable to NIST Standard reference materials.

Reference intervals

Concentrations are formulated at normal and abnormal levels expected in patient blood.

The expected clinical ranges of patient blood is referenced in Tietz, N.W. ed. 1995 Textbook of Clinical Chemistry, W.B. Saunders Co.

Users are encouraged to determine Expected Ranges in their own laboratory.²

Expected Ranges

The expected range for each analyte was determined at Nova Biomedical using deviations from the mean values that may be expected under differing laboratory conditions for analyzers from different manufacturers. The expected ranges for the Stat Profile pHox Ultra/Critical Care Xpress are determined by the manufacturer and are listed on the cartridge. Performance Check Solution results for glucose and creatinine will be Pass or Fail.

¹ NCCLS Document M29-T2.

² How to define and determine reference intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13.

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