REF 58395,58396

15°C 30°C CE

# Stat Profile Prime Plus® VET Calibrator Cartridge

Stat Profile Prime Plus® VET Kalibrierkassette, Φυσιγγιό βαθμονόμησης VET Stat Profile Prime Plus®, Cartucho de calibrador Stat Profile Prime Plus® VET, Cartouche de calibrage Stat Profile Prime Plus® VET, Cartuccia calibratore VET Stat Profile Prime Plus®, Cartucho de calibrador VET Stat Profile Prime Plus®, Stat Profile Prime Plus® VET kalibrátor patron, Stat Profile Prime Plus® VET בודק, בודק כוונת כיון, Stat Profile Prime Plus® VET キャリブレータカートリッジ, Stat Profile Prime Plus® VET 校准試剂盒

LOT 23340033

2025-05-13

CAL A	≥ 495 mL	CAL B	≥ 250 mL	CAL C	≥ 245 mL	CAL D	≥ 250 mL	CAL F	≥ 850 mL	SOLN R	≥ 500 mL
HCO <sub>3</sub>	11 mmol/L	pH	6.834 @ 37°C	HCO <sub>3</sub>	22 mmol/L	pH	7.354 @ 37°C	PO <sub>2</sub>	103.1 mmHg	KCl	2 mol/L
Na <sup>+</sup>	135.5 mmol/L	H <sup>+</sup>	146.55 nmol/L	Cl <sup>-</sup>	69.0 mmol/L	H <sup>+</sup>	44.26 nmol/L	PO <sub>2</sub>	13.71 kPa		
K <sup>+</sup>	4.0 mmol/L	Na <sup>+</sup>	77.0 mmol/L	Glu	200 mg/dL					!	WASTE CONTAINER W
Cl <sup>-</sup>	108.0 mmol/L	K <sup>+</sup>	10.0 mmol/L	Glu	11.10 mmol/L						
iCa	1.04 mmol/L	iCa	1.92 mmol/L	Lac	10.0 mmol/L						
iCa	4.16 mg/dL	iCa	7.68 mg/dL	Lac	89.0 mg/dL						
iMg	0.50 mmol/L	iMg	1.50 mmol/L	BUN	50 mg/dL						
iMg	1.22 mg/dL	iMg	3.63 mg/dL	BUN	17.9 mmol/L						
Glu	80 mg/dL	tHb	15.5 g/dL								
Glu	4.44 mmol/L	tHb	155 g/L								
Lac	2.0 mmol/L	tHb	9.6 mmol/L								
Lac	17.8 mg/dL										
BUN	10 mg/dL										
BUN	3.6 mmol/L										

## EN

**Description**  
The Cartridge contains the following aqueous solutions within individual bags: 5 Calibrators (Calibrator A, B, C, D and F), Reference Solution, and a flexible waste bag. Each bag includes a septum with a septum. The exposed bag fittings are arranged in a line along the rear of the cartridge box. The septa are pierced during the insertion of the cartridge into the analyzer. A charge kit consisting of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation.

**Intended Use**  
For the quantitative determination of Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, iCa, iMg, BUN (Urea), Glucose, Lactate, pH, PO<sub>2</sub>, PO<sub>3</sub>%, Hct, Hb, O<sub>2</sub>Hb, COHb, MetHb, HHb, and tHb in animal whole blood using the Stat Profile Prime Plus VET Analyzer.

**Methodology**  
Refer to Stat Profile Prime Plus VET Analyzer Instructions For Use Manual for Methodologies and Principles.

**Composition**  
Calibrator and solution volumes are listed in the table. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF. NCCLS DOCUMENT M29-T2).

**Warnings and Cautions**  
Refer to the Stat Profile Prime Plus VET Analyzer Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory reagents. Discard Cartridge in accordance with local guidelines.  
Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

**Storage**  
Store at 15 - 30°C. Expiration date is printed on the Calibrator Cartridge.

**Directions For Use**  
Refer to the Stat Profile Prime Plus VET Analyzer Instructions for Use for complete directions. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents.

**Traceability of Standards**  
Analytes traced to NIST Standard Reference Materials wherever possible.

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