Stat Profile® Prime Calibrator Cartridge CCS/CCS Comp

Kalibratorkassette, Καρέτα Βαθμονόμητής, Cartucho del calibrador, Cartouche étalon, Cartuccia calibratore, Calibrator Cartridge, Kalibratorkassett, Kalibratorkazetta, キャリブレータ カートリッジ, 측정기 카트리지, 校准液盒

LOT 23355035

2025-05-30

CAL A	≥ 500 mL	CAL B	≥ 250 mL	CAL F	≥ 720 mL	SOLN R	≥ 300 mL
pH	7.352	pH	6.840	PCO_2	52.2 mmHg	KCl	2 mol/L
PCO_2	24.9 mmHg	PO_2	167.0 mmHg	PO_2	97.2 mmHg		
Na^+	133.0 mmol/L	Na^+	68.5 mmol/L				
K^+	4.0 mmol/L	K^+	10.0 mmol/L				
Cl^-	106.0 mmol/L	Cl^-	46.0 mmol/L				
iCa	1.10 mmol/L	iCa	2.20 mmol/L				
Glu	80 mg/dL	Glu	200 mg/dL				
Lac	2.0 mmol/L	Lac	10.0 mmol/L				

(EN)

Product Description
The above table lists the concentrations of constituents for the calibrators and solutions contained in the Calibrator Cartridge. Contains the following aqueous solutions within individual bags housed in a cartridge box: 3 Calibrators (Calibrator A, B, and F), reference solution, and a flexible waste bag. Each bag includes a septum with a septa. The exposed bag flaps are arranged in a line along the rear of the cartridge box. These septa are pierced during the insertion of the cartridge into the analyzer.

Intended Use

The Stat Profile Prime Calibrator Cartridge CCS is intended for the calibration of pH, PCO_2 , PO_2 , Na^+ , K^+ , Cl^- , iCa, Glucose (Glu), and Lactate (Lac) using the Stat Profile Prime CCS Analyzer.

Methodology

Refer to Prime Analyzer Instructions For Use Manual for Test Methodologies and Principles.

Computation

Calibrator volumes are listed in the table above. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials, required for handling laboratory reagents. Discard Calibrator Cartridges and Waste Container after use.

Warnings and Cautions:

Intended for in vitro diagnostic use.

Refer to the Prime Analyzer Instructions for Use Manual for complete directions for use, installation, cleaning, recommendations for handling, and troubleshooting information. Methodology and Principles of the testing procedures. Follow standard practices required for handling laboratory reagents. Discard Calibrator Cartridges and Waste Container after use.

Storage

Store at 15–30°C. Expiration date is printed on the front of the Calibrator Cartridge. Directions For Use Refer to the Prime Analyzer Instructions for Use for complete directions for installation of the cartridge, and recommendations for Calibration Intervals. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents. Cartridge is ready for use.

Traceability of Standards

Analytes traced to NIST Standard Reference Materials.

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