Stat Profile® pHox Ultra/Critical Care Xpress Chemistry Controls & Performance Check Solution Auto-Cartridge

Chemiekontrollen und Leistungstest-Lösung Auto-QK-Kassette, Autómatu Káseta Xhímikón Elegktrón & Διαλύματος Elegkou Atipodósois,
Cartucho automático para Controles para química y Solución de verificación de comportamiento, Cartouche CQ de contrôles chimiques et de
solution de contrôle de linéarité, Cartuccia automatica controlli chimici e soluzione controllo prestazioni, Controlos de química & Cartucho
automático de Solução de Verificação de Desempenho, Autokassett med kemiska kontroller och prestandakontrollösning

LOT

24025032

2025-07-12

CONTROL	4	5	PERFORMANCE	CHECK	SOLN
----------------	----------	----------	--------------------	--------------	-------------

Expected Ranges, Erwartete Bereiche, Αναμενόμενα Εύρη, Rangos esperados, Fourchettes attendues, Intervalli previsti, Gamas previstas, Förväntade områden

		CONTROL 4	CONTROL 5
		min - \bar{x} - max	min - \bar{x} - max
Na ⁺	mmol/L	139.2 - 143.2 - 147.2	112.8 - 116.8 - 120.8
K ⁺	mmol/L	3.80 - 4.05 - 4.30	6.28 - 6.58 - 6.88
Cl ⁻	mmol/L	122.2 - 126.7 - 131.2	91.5 - 96.0 - 100.5
Ca ⁺⁺	mmol/L	1.03 - 1.11 - 1.19	1.52 - 1.62 - 1.76
Mg ⁺⁺	mmol/L	0.54 - 0.63 - 0.72	0.98 - 1.13 - 1.28
Glu	mg/dL	73 - 80 - 87	255 - 280 - 305
Glu	mmol/L	4.05 - 4.44 - 4.83	14.15 - 15.54 - 16.93
Lac	mmol/L	1.7 - 2.0 - 2.3	6.3 - 7.0 - 7.7
BUN	mg/dL	14 - 19 - 24	49 - 56 - 63
BUN	mmol/L	5.0 - 6.8 - 8.6	17.5 - 20.0 - 22.5
Urea	mg/dL	30.0 - 40.8 - 51.5	105.1 - 120.1 - 135.1
Urea	mmol/L	5.0 - 6.8 - 8.6	17.5 - 20.0 - 22.5

PERFORMANCE	CHECK	SOLN
PASS - No Action		
FAIL - Replace Membrane		
FEHLERFREI - Keine Maßnahme		
FEHLERHAFT - Membran austauschen		
ΕΠΙΤΥΧΙΑ - Καμία Ενέργεια		
ΑΠΟΤΥΧΙΑ - Αντικατοπτρίστε τη μεμβράνη		
PASÓ - Ninguna acción		
NO PASÓ - Sustituir membrana		
SUCCES - Aucune action		
ECHEC - Remplacer la membrane		
RIUSCITO - Nessuna azione		
NON RIUSCITO - Sostituisci membrana		
PASSOU - Nenhuma ação		
FALHOU - Substituir membrana		
OK - Ingen åtgärd		
FEL - Byt ut membran		

Product Description
Carton contains three bags. Two bags contain assayed quality control material formulated for monitoring the measurement of Na⁺, K⁺, Cl⁻, Ca⁺⁺, Mg⁺⁺, Glucose, Lactate, and BUN (Controls 4 and 5) and one bag contains an aqueous quality control material to verify glucose membrane performance (Performance Check Solution, Control 6).

Intended Use
For in vitro diagnostic use to monitor the performance of the Stat Profile pHox Ultra/Critical Care Xpress Analyzer and to insure optimum performance of glucose membranes.

For use in the laboratory, the solution should be analyzed before controls are run and after changing glucose membranes.

Methodology

Refer to Stat Profile pHox Ultra/Critical Care Xpress Analyzer Instruction for Use Manual.

Composition

Controls (Controls 4 and 5) contain known concentrations of Na⁺, K⁺, Cl⁻, Ca⁺⁺, Mg⁺⁺, Glucose, Lactate, and BUN at two levels. Performance Check Solution (Control 6) is a buffered salt solution containing Glucose (1000 mg/dL). Each bag contains a minimum volume of 90 mL.

Contaminants not of human origin, however, good laboratory practice should be followed during handling of these materials.

Warnings and Precautions

For in vitro diagnostic use. Do not freeze. Refer to the Stat Profile pHox Ultra/Critical Care Xpress Instructions for Use Manual for complete directions. Follow standard practices required during laboratory controls.

Storage

Store at 2-8°C. Do not freeze. Each cartridge has lot number and expiration date printed on the label.

Directions for Use

Refer to Analyzer Instructions for Use Manual for complete directions. Performance Check Solution and Quality Control cartridges are run and after changing glucose membranes. Verify that the lot number appearing on the expected ranges table is the same as the lot number printed on the cartridge.

Limitations

For use by Nova Biomedical Analyze Only. Performance characteristics have not been established for use on other manufacturer's analyzers.

Once installed, each Control Auto-Cartridge may be used for a maximum of 35 days from the date of initial activation and installation. The Analyzer will indicate that the cartridge is invalid at that time.

This product is warranted to perform as described in this product insert, and Nova Biomedical claims any implied warranty or merchantability or fitness for any other purpose, and in no event shall Nova Biomedical be liable for any consequential damages arising out of the aforesaid express warranties.

Traceability of Standards

Analyses are traced to NIST Standard reference materials.

Reference Intervals

Concentrations are formulated at normal and abnormal levels expected in patient blood. The expected clinical ranges of patient blood is referenced in Tietz, N.W. ed. 1986 Textbook of Clinical Biochemistry. Users may wish to determine Expected Ranges in their own laboratory.²

Expected Ranges

The expected range for each analyte was determined at Nova Biomedical using replicate determinations on the analyzer. The expected range is the range of values expected to deviate from the mean value, that may be experienced under different laboratory conditions for analyzers operating within specifications. Refer to Expected Ranges Table for Control 4 and 5 values.

Performance Check Solution results for glucose will be Pass or Fail.

¹ NCCLS Document M29-T2.

² Definition and determination of Reference intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13.

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