Stat Profile Prime Plus® Blood Gas, CO-Oximeter, Chemistry Controls Auto-Cartridge

Cartucho automático para controles de química, gases en sangre y cooxímetro Stat Profile Prime Plus®. Cartouche automatique de contrôles de gaz du sang/CO-oxymètre, chimie Stat Profile Prime Plus®, Stat Profile Prime Plus® Auto-Kassette für Blutgas-, CO-Oximeter- und Blutchemiekontrollen, Stat Profile Prime Plus® Αυτόματο φασίγγιο επιπέδου ελέγχου χημείας αερίων αίματος, CO-Οξόμετρο, Cartuccia per controlli automatici chimici per gas ematici/CO-ossimetria Stat Profile Prime Plus®, Cartucho automático de controles de química, de CO-oxímetero e de gas no sangue Stat Profile Prime Plus®, Stat Profile Prime Plus® vágáz/CO-oximéter, kémiai kontrollok automatikus patron, Stat Profile Prime Plus® מדידת בדיקות כימיה אוטומטית לזרם Co-Oximeter, Stat Profile Prime Plus® 血液ガス、CO オキシメーター、生化学検査用コントロール自動カートリッジ, Stat Profile Prime Plus® 혈액 가스, CO-산소 농도계, 화학 조제액 자동 카트리지, Stat Profile Prime Plus® 血气、一氧化碳-血氧仪、化学对照溶液自动试剂盒

LOT 24053066

CONTROL 1 2 3 4 5

2025-08-08

Table with 7 columns: Parameter, Unit, CONTROL 1 (min - X - max), CONTROL 2 (min - X - max), CONTROL 3 (min - X - max), CONTROL 4 (min - X - max), CONTROL 5 (min - X - max). Rows include pH, H+, PCO2, PCO2, PO2, PO2, SO2, Hct, Na+, K+, Cl-, iCa, iCa, iMg, iMg, Glu, Glu, Lac, Lac, HbF\*, iHb, iHb, iHb, O2Hb, COHb, MetHb, HbH, iBI1\*, iBI1\*, iBI1\*.

EN

Description
Aquous quality control material for monitoring the performance of pH, PCO2, PO2, SO2, hematocrit (Hct), fetal hemoglobin (HbF), total hemoglobin (HbT), total bilirubin (TBil), cyanmethemoglobin (COHb), methemoglobin (MetHb), and oxygen hemoglobin (O2Hb) in Levels 1, 2 and 3 as well as Na+, K+, Cl-, iCa, iMg, glucose, and lactate in Level 4 and 5. For use with Stat Profile Prime Plus Analyzers ONLY.

Intended Use
Intended for in vitro diagnostic use by healthcare professionals for monitoring the performance of the Stat Profile Prime Plus Analyzers.

Methology
Refer to Stat Profile Prime Plus Analyzer Instructions For Use Manual for Methodology and Principles.

Composition
Controls Levels 1, 2 and 3 are buffered bicarbonate solutions containing dye, salts and preservatives. Each level has a known pH and is equilibrated to a known O2, CO2, and N2 value. Controls Level 4 and 5 are buffered solutions containing known concentrations of Na+, K+, Cl-, iCa, iMg, glucose, lactate and preservatives. Each pouch contains a minimum of 100 mL. Controls contain no constituents of human origin, however good laboratory practices should be followed during handling of these materials.

Warnings and Cautions:
DO NOT FREEZE. Mix the cartridge by gently inverting for several seconds. DO NOT SHAKE CARTRIDGE. Refer to Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete information. Intended for in vitro diagnostic use. Follow standard practices for handling laboratory reagents.

Storage
Store at 2-8°C (37-48°F). DO NOT FREEZE.

Directions for use
Ensure controls are at room temperature prior to installation. Mix Cartridge well by gently inverting for 1 minute. Verify that the Lot Number on the Expected Ranges Table corresponds to the Lot Number on the cartridge. Refer to Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete directions.

Limitations
PO2 values vary inversely with temperature (approximately 1%/°C). Therefore, it is critical to follow the temperature guidelines described in "Directions for Use." The expected Range values are specific for instruments and controls manufactured by Nova Biomedical. Once installed, each Stat Profile Prime Plus Cartridge may be used for a maximum of 90 days from the initial installation date on the system at which time the system will indicate the cartridge is invalid. Each cartridge may be inserted and removed from the analyzer a maximum of 5 times.

Traceability of Standards
Total Hemoglobin (HbT) and Methemoglobin (MetHb) are traceable by using Cyanmethemoglobin method. Carboxyhemoglobin (COHb) and Oxygen hemoglobin (O2Hb) are traceable using Spectrophotometry. Analytes are traced to NIST Standard Reference Materials.

Reference Intervals
Concentrations are formulated at normal and abnormal expected values in patient blood. The expected clinical range of these values as patient blood is referenced in Tietz, WJ ed. 1988 Textbook of Clinical Chemistry, W.B. Saunders Co. Users may wish to determine Mean Values and Expected Ranges in their own laboratory.

Expected Ranges
The expected range for each parameter was determined at Nova Biomedical using replicate determinations on Nova analyzers. The expected range indicates the maximum deviations from the Mean Value that may be expected under differing laboratory conditions for instruments operating within specifications. Refer to Expected Ranges Table.

\*Not available in the USA or for Point-of-Care/Non-Patient Testing use.

INCLOS Document M29-12
Zřizov a Definice a Deternace Referenční Intervalů in the clinical laboratory, approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13.

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