Stat Profile® Prime Calibrator Cartridge CCS/CCS Comp

Kalibratorkassette, Κασέτα Βαθμονομητή, Cartucho del calibrador, Cartouche étalon, Cartuccia calibratore, Calibrator Cartridge, Kalibratorkassett, Kalibratorkazetta, キャリブレータ カートリッジ, 측정기 카트리지, 校准液盒

LOT 24046062

2025-08-01

CAL A	≥ 500 mL	CAL B	≥ 250 mL	CAL F	≥ 720 mL	SOLN R	≥ 300 mL
pH	7.359	pH	6.840	PCO_2	52.2 mmHg	KCl	2 mol/L
PCO_2	24.4 mmHg	PO_2	167.0 mmHg	PO_2	95.7 mmHg		
Na^+	133.0 mmol/L	Na^+	68.5 mmol/L				
K^+	4.0 mmol/L	K^+	10.0 mmol/L				
Cl^-	106.0 mmol/L	Cl^-	46.0 mmol/L				
iCa	1.10 mmol/L	iCa	2.20 mmol/L				
Glu	80 mg/dL	Glu	200 mg/dL				
Lac	2.0 mmol/L	Lac	10.0 mmol/L				

(EN)

Product Description

The above table lists the concentrations of constituents for the calibrators and solutions contained in the Stat Profile Prime Calibrator Cartridge CCS. The calibrators and solutions included in the table are kept in a cardboard box: 3 Calibrators (Calibrator A, B and F), reference solution, and a flexible waste bag. Each bag includes a septum with a septa. The exposed bag fittings are arranged in a line along the rear of the cartridge box. These septa are pierced during the insertion of the cartridge into the analyzer.

Intended Use

The Stat Profile Prime Calibrator Cartridge CCS is intended for the calibration of pH, PCO_2 , PO_2 , Hct , Na^+ , K^+ , iCa , Glucose (Glu), and Lactate (Lac) using the Stat Profile Prime CCS analyser.

Methodology

Refer to Prime Analyzer Instructions For Use Manual for Test Methodologies and Principles.

Composition

Calibrator volumes are listed in the table above. Contains no constituents of human origin, however good laboratory practice should be followed during handling of these materials. (REF. NCCLS DOCUMENT M29-T2).

Warnings and Cautions:

Intended for in vitro diagnostic use.

Refer to the Prime Analyzer Instructions For Use Manual for complete directions for use, installation instructions, and recommendations for calibration intervals. Calibrator and solution components are supplied in individual sealed bags. Do not reuse after opening for handling laboratory reagents. Discard Calibrator Cartridges and Waste products in accordance with local guide lines. Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

Storage

Store at 15-30°C. Expiration date is printed on the front of the Calibrator Cartridge.

Directions For Use Refer to the Prime Analyzer Instructions For Use for complete directions for installation of the cartridge, and recommendations for Calibration Intervals. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents. Cartridge is ready for use.

Traceability of Standards

Analytes traced to NIST Standard Reference Materials.

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