REF 48927, 48928

IVD 2°C-8°C CE

# Stat Profile® pHox Ultra/Critical Care Xpress Chemistry Controls & Performance Check Solution Auto-Cartridge With Creatinine

Chemiekontrollen und Leistungstest-Lösung Auto-QK-Kassette mit Creatinin, Αυτόματη Κασέτα Χημικών Ελεγκτών & Διαλύματος Ελέγχου Απόδοσης Με Κρεατινίνη, Cartucho automático para Controles para química y Solución de verificación de comportamiento con creatinina, Cartouche CQ de contrôles chimiques et de solution de contrôle de linéarité avec créatinine, Cartuccia automatica controlli chimici e soluzione controllo prestazioni con creatinina, Controles de química & Cartucho automático de Solução de Desempenho com creatinina, Autokassett med kemiska kontroll och prestandakontrollösning med kreatinin

LOT 24057037 2025-07-17

CONTROL 4 5 PERFORMANCE CHECK SOLN

## Expected Ranges, Erwartete Bereiche, Αναμενόμενα Εύρη, Rangos esperados, Fourchettes attendues, Intervalli previsti, Gamas previstas, Förväntade områden

		CONTROL 4			CONTROL 5		
		min	$\bar{x}$	max	min	$\bar{x}$	max
Na <sup>+</sup>	mmol/L	136.0	140.0	144.0	110.5	114.5	118.5
K <sup>+</sup>	mmol/L	3.72	3.97	4.22	6.10	6.40	6.70
Cl <sup>-</sup>	mmol/L	121.3	125.8	130.3	91.5	96.0	100.5
Ca <sup>++</sup>	mmol/L	1.03	1.11	1.19	1.51	1.63	1.75
Mg <sup>++</sup>	mmol/L	0.55	0.64	0.73	1.00	1.15	1.30
Glu	mg/dL	73	80	87	255	280	305
Glu	mmol/L	4.05	4.44	4.83	14.15	15.54	16.93
Lac	mmol/L	1.7	2.0	2.3	6.3	7.0	7.7
BUN	mg/dL	14	19	24	49	56	63
BUN	mmol/L	5.0	6.8	8.6	17.5	20.0	22.5
Urea	mg/dL	30.0	40.8	51.5	105.1	120.1	135.1
Urea	mmol/L	5.0	6.8	8.6	17.5	20.0	22.5
Creatinine	mg/dL	0.7	1.0	1.3	5.2	6.2	7.2
Creatinine	μmol/L	61.9	88.5	115.0	460.0	548.5	637.0

	PERFORMANCE CHECK		SOLN
	PASS - No Action	FAIL - Replace Membrane	
GLUCOSE	FEHLERFREI - Keine Maßnahme	FEHLERHAFT - Membran austauschen	CREATININE
	ΕΠΙΤΥΧΙΑ - Καμία Ενέργεια	ΕΠΙΤΥΧΙΑ - Αντικαταστήστε τη μεμβράνη	
	PASO - Ninguna acción	NO PASO - Sustituir membrana	
	SUCCES - Aucune action	ECHEC - Remplacer la membrane	
	RIUSCITO - Nessuna azione	NON RIUSCITO - Sostituisce membrana	
	PASSOU - Nenhuma ação	FALHOU - Substituir membrana	
	OK - Ingen åtgärd	FEL - Byt ut membran	
	FEHLERFREI - Keine Maßnahme	FEHLERHAFT - Membran austauschen	
	ΕΠΙΤΥΧΙΑ - Καμία Ενέργεια	ΕΠΙΤΥΧΙΑ - Αντικαταστήστε τη μεμβράνη	
	PASO - Ninguna acción	NO PASO - Sustituir membrana	
SUCCES - Aucune action	ECHEC - Remplacer la membrane		
RIUSCITO - Nessuna azione	NON RIUSCITO - Sostituisce membrana		
PASSOU - Nenhuma ação	FALHOU - Substituir membrana		
OK - Ingen åtgärd	FEL - Byt ut membran		

Product Description  
Carton contains three bags. Two bags contain assayed quality control material formulated for monitoring the measurement of Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, Ca<sup>++</sup>, Mg<sup>++</sup>, Glucose, Lactate, BUN and Creatinine (Control 4 and 5). One bag contains an aqueous quality control material to verify glucose and creatinine membrane performance (Performance Check Solution).

Intended Use  
For in vitro diagnostic use to monitor the performance of the Stat Profile pHox Ultra/Critical Care Xpress Analyzer and to measure optimum performance of glucose and creatinine membranes. Performance check solution should be analyzed before controls are run and after changing glucose or creatinine membranes.

Methodology  
Refer to Stat Profile pHox Ultra/Critical Care Xpress Analyzer Instructions for Use Manual.

Composition  
Controls (4 and 5) contain known concentrations of Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, Ca<sup>++</sup>, Mg<sup>++</sup>, Glucose, Lactate, BUN, and Creatinine at two levels. Performance Check Solution is a buffered salt solution containing known levels of Glucose, Creatinine, and Creatine. Each bag contains a minimum volume of 100 mL.

Warnings and Precautions  
For in vitro diagnostic use. Do not freeze. Refer to the Stat Profile pHox Ultra/Critical Care Xpress Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory controls.

Storage  
Store at 2-8°C. Do not freeze. Each cartridge has lot number and expiration date printed on the label.

Directions for Use  
Activation of Cartridge  
Each level of control must be charged with the creatinine syringe prior to installation of the Controls Auto-Cartridge. The three syringes are labeled and color coded to correspond to the firments on the carton for each level of control.

Limitations  
Use on Nova Biomedical Analyzers Only. Once installed, the Auto-Cartridge may be used for a maximum of 21 days from the date of initial activation and installation. The Analyzer will indicate that the cartridge is invalid at that time.

Traceability of Controls/Solution  
Analyses are traced to NIST Standard reference materials.

Reference Intervals  
Concentrations are formulated at normal and abnormal levels expected in patient blood.

Expected Ranges  
The expected range for each analyte was determined at Nova Biomedical using replicate determinations on Nova Analyzers. The expected range indicates the maximum deviations from the mean value that may be expected under differing laboratory conditions for analyzers operating within specifications. Refer to Expected Ranges Table for Control 4 and 5 values. Performance Check Solution results for glucose and creatinine will be Pass or Fail.

How to define and determine reference intervals in the clinical laboratory: approved guideline - second edition, NCLLS C28-A2, Volume 20, Number 13.

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**EC REP** Nova Biomedical GmbH
Waltham, MA 02454 U.S.A.