# Stat Profile Prime Plus® Calibrator Cartridge

Stat Profile Prime Plus® Kalibrerkassette, Φυσιγγο μέσου βαθμονόμησης με κρεατηνή Stat Profile Prime Plus®, Cartucho de calibrador Stat Profile Prime Plus®, Cartouche de calibrage Stat Profile Prime Plus®, Cartuccia calibratore Stat Profile Prime Plus®, Cartucho de calibrador Stat Profile Prime Plus®, Stat Profile Prime Plus® kalibrátor patron, 专利权 Stat Profile Prime Plus®, Stat Profile Prime Plus® キャリブレータ カートリッジ, Stat Profile Prime Plus® 보정기 카트리지, Stat Profile Prime Plus® 校准试剂盒

LOT 24059112

2025-08-16

CAL A	CAL B	CAL C	CAL D	CAL F	SOLN R
HCO <sub>3</sub> <sup>-</sup> 11 mmol/L	pH 6.834 @ 37°C	HCO <sub>3</sub> <sup>-</sup> 22 mmol/L	pH 7.354 @ 37°C	PO <sub>2</sub> 104.4 mmHg	KCl 2 mol/L
Na <sup>+</sup> 135.5 mmol/L	H <sup>+</sup> 146.55 nmol/L	Cl <sup>-</sup> 69.0 mmol/L	H <sup>+</sup> 44.26 nmol/L	PO <sub>2</sub> 13.89 kPa	
K <sup>+</sup> 4.0 mmol/L	Na <sup>+</sup> 77.0 mmol/L	Glu 200 mg/dL			
Cl <sup>-</sup> 108.0 mmol/L	K <sup>+</sup> 10.0 mmol/L	Glu 11.10 mmol/L			
iCa 1.04 mmol/L	iCa 1.92 mmol/L	Lac 10.0 mmol/L			
iCa 4.16 mg/dL	iCa 7.68 mg/dL	Lac 89.0 mg/dL			
iMg 0.50 mmol/L	iMg 1.50 mmol/L	BUN 50 mg/dL			
iMg 1.22 mg/dL	iMg 3.63 mg/dL	BUN 17.9 mmol/L			
Glu 80 mg/dL	tHb 15.5 g/dL				
Glu 4.44 mmol/L	tHb 155 g/L				
Lac 2.0 mmol/L	tHb 9.6 mmol/L				
Lac 17.8 mg/dL					
BUN 10 mg/dL					
BUN 3.6 mmol/L					

! WASTE CONTAINER W

## EN

Description  
The Cartridge contains the following aqueous solutions within individual bags: 5 Calibrators (Calibrator A, B, C, D, and F SOLN), Reference Solution, and a flexible waste bag. Each bag includes a fitting with a septum. The exposed bag fittings are arranged in a line along the rear of the cartridge box. The septa are pierced during the insertion of the cartridge into the analyzer.

Intended Use  
For the quantitative determination of Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, iCa, iMg, BUN (Urea), Glucose, Lactate, pH, PCO<sub>2</sub>, PO<sub>2</sub>, SO<sub>4</sub><sup>2-</sup>, Hct/Hb, Hbf<sup>-</sup>, O<sub>2</sub>Hb, COHb, MetHb, HbB, and tBil<sup>+</sup> in human blood using the Stat Profile Prime Plus Analyzer.

Not available in the USA or for Point-of-Care/Near-Patient Testing use.

Methodology  
Refer to Stat Profile Prime Plus Analyzer Instructions For Use Manual for Methodologies and Principles.

Composition  
The calibrators are aqueous salt solutions containing analytes at concentrations listed in the table above. Calibrator B also contains dye. Contain no constituents of human origin, however, good laboratory practice should be followed during the handling of these materials. (REF. NCCLS DOCUMENT M29-T2).

Warnings and Cautions  
Intended for in vitro diagnostic use.  
Refer to the Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory reagents. Discard Cartridge in accordance with local guidelines.

Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

Storage  
Store at 15 - 30°C. Expiration date is printed on the Calibrator Cartridge.

Directions For Use  
Refer to the Stat Profile Prime Plus Analyzer Instructions for Use for complete directions. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents.

Traceability of Standards  
Analytes traced to NIST Standard Reference Materials wherever possible.

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