REF 57825, 57826, 57827, 57828, 57829, 65141

Rx Only IVD 15°C 30°C CE

# Stat Profile Prime Plus® Calibrator Cartridge

Stat Profile Prime Plus® Kaliibratorkassette, Φυσιγγιο μέσου βαθμονόμησης με κραστινή Stat Profile Prime Plus®, Cartucho de calibrador Stat Profile Prime Plus®, Cartouche de calibrateurs Stat Profile Prime Plus®, Cartuccia calibratore Stat Profile Prime Plus®, Cartucho de calibrador Stat Profile Plus®, Stat Profile Prime Plus® kalibrátor patron, גל' ציינית Stat Profile Prime Plus®, Stat Profile Prime Plus® Stat Profile Prime Plus® カリブレーターカートリッジ, Stat Profile Prime Plus® 보정기 카트리지, Stat Profile Prime Plus® 校准试剂盒

LOT 24065074

2025-07-18

CAL A	CAL B	CAL C	CAL D	CAL F	SOLN R
HCO <sub>3</sub> <sup>-</sup> 11 mmol/L	pH 6.834 @ 37°C	HCO <sub>3</sub> <sup>-</sup> 22 mmol/L	pH 7.354 @ 37°C	Po <sub>2</sub> 104.1 mmHg	KCl 2 mol/L
Na <sup>+</sup> 135.5 mmol/L	H <sup>+</sup> 146.55 nmol/L	Cl <sup>-</sup> 69.0 mmol/L	H <sup>+</sup> 44.26 nmol/L	Po <sub>2</sub> 13.85 kPa	WASTE CONTAINER W
K <sup>+</sup> 4.0 mmol/L	Na <sup>+</sup> 77.0 mmol/L	Glu 200 mg/dL			
Cl <sup>-</sup> 108.0 mmol/L	K <sup>+</sup> 10.0 mmol/L	Glu 11.10 mmol/L			
iCa 1.04 mmol/L	iCa 1.92 mmol/L	Lac 10.0 mmol/L			
iCa 4.16 mg/dL	iCa 7.68 mg/dL	Lac 89.0 mg/dL			
iMg 0.50 mmol/L	iMg 1.50 mmol/L	BUN 50 mg/dL			
iMg 1.22 mg/dL	iMg 3.63 mg/dL	BUN 17.9 mmol/L			
Glu 80 mg/dL	tHb 15.4 g/dL				
Glu 4.44 mmol/L	tHb 154 g/L				
Lac 2.0 mmol/L	tHb 9.6 mmol/L				
Lac 17.8 mg/dL					
BUN 10 mg/dL					
BUN 3.6 mmol/L					

## EN

## Description

The Cartridge contains the following aqueous solutions within individual bags: 5 Calibrators (Calibrator A, B, C, D, and F SOLN), Reference Solution, and a flexible waste bag. Each bag contains a filament with a septum. The exposed bag filaments are arranged in a line along the top edge of the cartridge box. The septa are pierced during the insertion of the cartridge into the analyzer.

## Intended Use

For the quantitative determination of Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, iCa, iMg, BUN (Urea), Glucose, Lactate, pH, Po<sub>2</sub>, Po<sub>2</sub> %, Hct, tHb, Hbf, O<sub>2</sub>Hb, COHb, MetHb, HHb, and iMg<sup>2+</sup> in human blood using the Stat Profile Prime Plus Analyzer.

\*Not available in the USA or for Point-of-Care/Near-Patient Testing use.

## Methodology

Refer to the Stat Profile Prime Plus Analyzer Instructions For Use Manual for Methodologies and Principles.

## Composition

The calibrators are aqueous salt solutions containing analytes at concentrations listed in the table above. Calibrator B also contains dye. Contain no constituents of human origin, however, standard laboratory practices should be followed during the handling of these materials. (REF. NCCLS DOCUMENT M29-T2).

## Warnings and Cautions

Intended for In vitro diagnostic use.

Refer to the Stat Profile Prime Plus Analyzer Instructions For Use Manual for complete directions. Follow standard practices required for handling laboratory reagents. Discard Cartridge in accordance with local guidelines.

Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

## Storage

Store at 15 - 30°C. Expiration date is printed on the Calibrator Cartridge.

## Directions For Use

Refer to the Stat Profile Prime Plus Analyzer Instructions For Use for complete directions. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents.

## Traceability of Standards

Analytes traced to NIST Standard Reference Materials wherever possible.

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EC REP Nova Biomedical GmbH Hessenring 13A, Geb. G 64546 Mörfelden-Walldorf Germany
Waltham, MA 02454 U.S.A.