Stat Profile Prime Plus® VET Blood Gas, CO-Oximeter, Chemistry Controls Auto-Cartridge with Creatinine

Cartucho automático con creatinina para controles de química, gases en sangre y cooxímetro Stat Profile Prime Plus® VET, Cartouche automatique de contrôles de gaz du sang/CO-oxymètre, chimie Stat Profile Prime Plus® VET avec créatine, Stat Profile Prime Plus® VET Auto-Kassette mit Kreatinin für Blutgas, CO-Oximeter- und Blutchemiekontrollen, Αυτόματο φορτίσιμο με κρεατινίνη επιπέδου ελέγχου χημείας αερίων αίματος, CO-Οξυμέτρο, CO-Οξυμέτρο Stat Profile Prime Plus® VET, Cartuccia con creatinina per controlli automatici chimici per gas ematici/CO-ossimetria Stat Profile Prime Plus® VET, Cartucho automático de controles de química, de CO-oxímetro e de gás no sangue Stat Profile Prime Plus® VET com creatinina, Stat Profile Prime Plus® VET vőrgáz, CO-oximéter, kémiai kontrollok automatikus patron kreatininnel, (מחשבונית בקרת אוויר ומסתה עם קריאטינין) Stat Profile Prime Plus® VET Chemistry Controls Auto-Cartridge, Stat Profile Prime Plus® VET 血液ガス, CO オキシメーター, 生化学検査用コントロール自動カートリッジ(クレアチニン), 크레아티닌 사용 Stat Profile Prime Plus® VET 혈액 가스, CO-산소 농도계, 화학 조절제 자동 카트리지, Stat Profile Prime Plus® VET 血气、一氧化碳-血氧仪、化学对照溶液自动试剂盒 (含肌酐)

LOT 24064040

CONTROL 1 2 3 4 5

2025-08-22

Table with 7 columns: Parameter, CONTROL 1, CONTROL 2, CONTROL 3, CONTROL 4, CONTROL 5. Rows include pH, H+, PCO2, PCO2, PO2, PO2, SO2, Hct, Na+, K+, Cl-, iCa, iCa, iMg, iMg, Glu, Glu, Lac, Lac, BUN, BUN, Urea, Urea, Creatinine, Creatinine, Creatinine, HbF, tHb, tHb, tHb, O2Hb, COHb, MetHb, HHb, tBil, tBil.

EN

Product Description
Aquinos control material for monitoring the performance of pH, PCO2, PO2, SO2, hematoct (hct), fetal hemoglobin (HbF), total hemoglobin (tHb), carboxyhemoglobin (COHb), methemoglobin (MetHb), and deoxyhemoglobin (iHb) in Levels 1, 2, 3 and 5 as well as Na+, K+, Cl-, Ca, Mg, glucose, lactate, BUN (urea), and Creatinine in Levels 4 and 5. A charge kit consisting of 12 syringes filled with creatinine solution is included for cartridge activation prior to installation. For use with Stat Profile Prime Plus VET Analyzers ONLY.

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Reference Intervals
Concentrations are formulated at normal and abnormal expected values in patient blood. The expected clinical range of these values in patient blood is referenced in Tietz, NW ed., 1986 Textbook of Clinical Chemistry, W.B. Saunders Co. Users may wish to determine Mean Values and Expected Ranges in their own laboratory?

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200 Prospect Street Waltham, MA 02454 U.S.A.