Stat Profile® Prime Calibrator Cartridge ABG

Kalibratorkassette, Κασέτα Βαθμονομής, Cartucho del calibrador, Cartouche étalon, Cartuccia calibratore, Calibrator Cartridge, Kalibratorkassett

LOT 24066046

2025-08-27

CAL A	$\geq 500 \text{ mL}$	CAL B	$\geq 250 \text{ mL}$	CAL F	$\geq 720 \text{ mL}$	SOLN R	$\geq 300 \text{ mL}$
pH	7.356	pH	6.840	PCO_2	51.0 mmHg	KCl	2 mol/L
PCO_2	24.5 mmHg	PO_2	167.0 mmHg	PO_2	98.0 mmHg	!	WASTE CONTAINER W

Product Description

The above table lists the concentrations of constituents for the calibrators and solutions contained in the Calibrator Cartridge. Contains the following aqueous solutions within individual flexible bags housed in a cardboard box: 3 Calibrators (Calibrator A, B, and F), Reference solution, and a flexible waste bag. Each bag includes a fitting with a septa. The exposed bag fittings are arranged in a line along the rear of the cartridge box. These septa are pierced during the insertion of the cartridge into the analyzer.

Intended Use

For the quantitative determination of pH, PCO_2 , and PO_2 in human blood using the Nova Biomedical Prime ABG Analyzer.

Methodology

Refer to the Prime Analyzer Instructions For Use Manual for Test Methodologies and Principles.

Composition

Calibrator volumes are listed in the table above. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF. NCCLS DOCUMENT M29-T2).

Warnings and Cautions:

Intended for *in vitro* diagnostic use.
Refer to the Prime Analyzer Instructions for Use Manual for complete directions for use, including instructions, recommendations for calibration intervals, Troubleshooting information, Methodology and Principles of the testing procedures. Follow standard practices required for handling laboratory reagents. Discard Calibrator Cartridges and Waste products in accordance with local guide lines. Waste Bag Contains Preservative (eyes, skin, and respiratory irritant).

Storage

Store at 15-30°C. Expiration date is printed on the front of the Calibrator Cartridge. Direct from the manufacturer. Instructions for Use for complete directions for installation of the cartridge, and recommendations for Calibration Intervals. Calibrator Cartridge must be inserted gently several times to thoroughly mix contents. Cartridge is ready for use.

Traceability of Standards

Analytes traced to NIST Standard Reference Materials.

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