REF 57831, 57832, 57833, 57834, 57835, 65142

Rx Only IVD 15°C 30°C CE

Stat Profile Prime Plus® Calibrator Cartridge With Creatinine

Stat Profile Prime Plus® Καλιβριεστική με Κρεατίνη, Stat Profile Prime Plus® Φυσιόγιο μέσω βραχυμόρφησης με κρεατίνη, Cartucho de calibrador con creatinina Stat Profile Prime Plus®, Cartouche de calibrateurs Stat Profile Prime Plus® avec créatinine, Cartuccia calibratore con creatinina Stat Profile Prime Plus®, Cartucho de calibrador Stat Profile Prime Plus® com creatinina, Stat Profile Prime Plus® kalibrátor patron kreatininnel, סט פרופיל פריים פלוס עם קריאטיν, Stat Profile Prime Plus® סט פרופיל פריים פלוס קריאטיν (קריאטיν), 크레아티닌, 크레아티닌 사용 Stat Profile Prime Plus® 보정기 카트리지, Stat Profile Prime Plus® 含肌酐的校准试剂盒

LOT 24092058

2025-09-19

CAL A		CAL B		CAL C		CAL D		CAL F		SOLN R	
HCO ₃	11 mmol/L	pH	6.834 @ 37°C	HCO ₃	22 mmol/L	pH	7.354 @ 37°C	PO ₂	99.1 mmHg	KCl	2 mol/L
Na ⁺	135.5 mmol/L	H ⁺	146.55 nmol/L	Cl ⁻	69.0 mmol/L	H ⁺	44.26 nmol/L	PO ₂	13.18 kPa	 WASTE CONTAINER W	
K ⁺	4.0 mmol/L	Na ⁺	77.0 mmol/L	Glu	200 mg/dL	Glu	11.10 mmol/L				
Cl ⁻	108.0 mmol/L	K ⁺	10.0 mmol/L	Glu	11.10 mmol/L	Lac	10.0 mmol/L				
iCa	1.04 mmol/L	iCa	1.92 mmol/L	Lac	10.0 mmol/L	Lac	89.0 mg/dL				
iCa	4.16 mg/dL	iCa	7.68 mg/dL	BUN	50 mg/dL	BUN	17.9 mmol/L				
iMg	0.50 mmol/L	iMg	1.50 mmol/L	BUN	17.9 mmol/L	Creat	5.0 mg/dL				
iMg	1.22 mg/dL	iMg	3.63 mg/dL	Creat	5.0 mg/dL	Creat	0.44 mmol/L				
Glu	80 mg/dL	tHb	15.4 g/dL	Creat	0.44 mmol/L	Creat	442 μmol/L				
Glu	4.44 mmol/L	tHb	154 g/L								
Lac	2.0 mmol/L	tHb	9.6 mmol/L								
Lac	17.8 mg/dL										
BUN	10 mg/dL										
BUN	3.6 mmol/L										
Creat	1.0 mg/dL										
Creat	0.09 mmol/L										
Creat	88.4 μmol/L										

EN
Description
The cartridge contains the following aqueous solutions within individual bags: 5 Calibrators (Calibrator A, B, C, D, and F), Reference Solution, and a flexible waste bag. Each bag includes a fitter with a septum. The exposed bag fitters are arranged in a line along the rear of the cartridge box. The septa are pierced during the insertion of the cartridge into the analyzer. A charge kit consisting of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation.

Intended Use
For the quantitative determination of Na⁺, K⁺, Cl⁻, iCa, iMg, BUN (Urea), Glucose, Lactate, Creatinine, pH, PCO₂, PO₂, SO₂, Hct, Hb, HbF, O₂Hb, COHb, Methb, HbH, and tHb⁺ in human blood using the Stat Profile Prime Plus Analyzer.

Methodology
Refer to Stat Profile Prime Plus Analyzer Instructions For Use Manual for Methodologies and Principles.

Composition
The calibrators are aqueous salt solutions containing analytes at concentrations listed in the table above. Calibrator B also contains dye. Contain no constituents of human origin, however, good laboratory practice should be followed during the handling of these materials. (REF: NCCLS DOCUMENT M29-T2).

Warnings and Cautions
Intended for in vitro diagnostic use. Refer to the Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory reagents. Discard Cartridge in accordance with local guidelines. Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

Storage
Store at 15 - 30°C. Expiration date is printed on the Calibrator Cartridge.

Directions For Use
Activation of Cartridge
Cal A and Cal C must be charged with creatinine syringe prior to installation of the Calibrator Cartridge. The two syringes are labeled and color coded to correspond to the fittings on the carton for Cal A and Cal C.

- Hold the syringe with tip side down and remove protective cap.
- Attach one of the enclosed needle assemblies to the syringe. Remove the protective cover from the needle.
- Match the color and label of the syringe to the appropriate fitter and insert needle.
- Slowly depress syringe plunger until the contents are dispensed.
- DO NOT PULL BACK ON THE PLUNGER TO FLUSH CONTENTS OF SYRINGE.
- Remove needle/syringe assembly from fitter and discard in an appropriate sharps container.
- Repeat Steps 1-5 for the next Cal.
- Mix Cartridge well by gently shaking for 1 minute.

Refer to the Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete directions. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents. Cartridge is ready for use.

*Not available in the USA or for Point-of-Care/Near-Patient Testing use.

Traceability of Standards
Analytes traced to NIST Standard Reference Materials wherever possible.

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EC REP Nova Biomedical GmbH Hessenring 13A, Geb. G 64546 Mörfelden-Walldorf Germany
Waltham, MA 02454 U.S.A.