REF 57825, 57826, 57827, 57828, 57829, 65141

Rx Only IVD 15°C - 30°C CE

# Stat Profile Prime Plus® Calibrator Cartridge

Stat Profile Prime Plus® Καλιβριζοκασετίδα, Φυσιόγιο μέσου βαθμονόμησης με κρεατινίνη Stat Profile Prime Plus®, Cartucho de calibrador Stat Profile Prime Plus®, Cartouche de calibrateurs Stat Profile Prime Plus®, Cartuccia calibratore Stat Profile Prime Plus®, Cartucho de calibrador Stat Profile Prime Plus®, Stat Profile Prime Plus® kalibrátor patron, כלי מחסנית Stat Profile Prime Plus®, Stat Profile Prime Plus® キャリブレーター カートリッジ, Stat Profile Prime Plus® 보정기 카트리지, Stat Profile Prime Plus® 校正试剂盒

LOT 24106045

2025-10-01

| CAL A            |              | CAL B           |               | CAL C            |              | CAL D          |              | CAL F           |            | SOLN R  |         |
|------------------|--------------|-----------------|---------------|------------------|--------------|----------------|--------------|-----------------|------------|---|---------|
| HCO <sub>3</sub> | 11 mmol/L    | pH              | 6.834 @ 37°C  | HCO <sub>3</sub> | 22 mmol/L    | pH             | 7.354 @ 37°C | PO <sub>2</sub> | 101.5 mmHg | KCl   | 2 mol/L |
| Na <sup>+</sup>  | 135.5 mmol/L | H <sup>+</sup>  | 146.55 nmol/L | Cl <sup>-</sup>  | 69.0 mmol/L  | H <sup>+</sup> | 44.26 nmol/L | PO <sub>2</sub> | 13.50 kPa  | <div style="border: 1px solid black; padding: 5px; text-align: center;"> <b>WASTE CONTAINER</b> <b>W</b> </div> |         |
| K <sup>+</sup>   | 4.0 mmol/L   | Na <sup>+</sup> | 77.0 mmol/L   | Glu              | 200 mg/dL    |                |              |                 |            |   |         |
| Cl <sup>-</sup>  | 108.0 mmol/L | K <sup>+</sup>  | 10.0 mmol/L   | Glu              | 11.10 mmol/L |                |              |                 |            |   |         |
| iCa              | 1.04 mmol/L  | iCa             | 1.92 mmol/L   | Lac              | 10.0 mmol/L  |                |              |                 |            |   |         |
| iCa              | 4.16 mg/dL   | iCa             | 7.68 mg/dL    | Lac              | 89.0 mg/dL   |                |              |                 |            |   |         |
| iMg              | 0.50 mmol/L  | iMg             | 1.50 mmol/L   | BUN              | 50 mg/dL     |                |              |                 |            |   |         |
| iMg              | 1.22 mg/dL   | iMg             | 3.63 mg/dL    | BUN              | 17.9 mmol/L  |                |              |                 |            |   |         |
| Glu              | 80 mg/dL     | tHb             | 15.4 g/dL     |                  |              |                |              |                 |            |   |         |
| Glu              | 4.44 mmol/L  | tHb             | 154 g/L       |                  |              |                |              |                 |            |   |         |
| Lac              | 2.0 mmol/L   | tHb             | 9.6 mmol/L    |                  |              |                |              |                 |            |   |         |
| Lac              | 17.8 mg/dL   |                 |               |                  |              |                |              |                 |            |   |         |
| BUN              | 10 mg/dL     |                 |               |                  |              |                |              |                 |            |   |         |
| BUN              | 3.6 mmol/L   |                 |               |                  |              |                |              |                 |            |   |         |

**EN**

**Description**  
The Cartridge contains the following aqueous solutions within individual bags: 5 Calibrators (Calibrator A, B, C, D, and F SOLN), Reference Solution, and a flexible waste bag. Each bag includes a filament with a septum. The exposed bag filaments are arranged in a line along the rear of the cartridge box. The septa are pierced during the insertion of the cartridge into the analyzer.

**Intended Use**  
For the quantitative determination of Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, iCa, iMg, BUN (Urea), Glucose, Lactate, pH, PCO<sub>2</sub>, PO<sub>2</sub>, SO<sub>2</sub>, Hct, Hb, HbF, O<sub>2</sub>Hb, COHb, MetHb, Hb, and tBil<sup>†</sup> in human blood using the Stat Profile Prime Plus Analyzer.

\*Not available in the USA or for Point-of-Care/Near-Patient Testing use.

**Methodology**  
Refer to Stat Profile Prime Plus Analyzer Instructions For Use Manual for Methodologies and Principles.

**Composition**  
The calibrators are aqueous salt solutions containing analytes at concentrations listed in the table above. Calibrator B also contains dye. Contain no constituents of human origin, however, good laboratory practice should be followed during the handling of these materials. (REF. NCCLS DOCUMENT M29-T2).

**Warnings and Cautions**  
Intended for *in vitro* diagnostic use.  
Refer to the Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory reagents. Discard Cartridge in accordance with local guidelines.  
Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

**Storage**  
Store at 15 - 30°C. Expiration date is printed on the Calibrator Cartridge.

**Directions For Use**  
Refer to the Stat Profile Prime Plus Analyzer Instructions for Use for complete directions. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents.

**Traceability of Standards**  
Analytes traced to NIST Standard Reference Materials wherever possible.

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