# Stat Profile Prime Plus® VET Calibrator Cartridge With Creatinine

Stat Profile Prime Plus® VET Kalibrierkassette mit Kreatinin, Φυσίγγιο μέσου βαθμονόμησης VET με κρεατινίνη Stat Profile Prime Plus®, Cartucho de calibrador con creatinina Stat Profile Prime Plus® VET, Cartouche de calibreurs avec creatinine Stat Profile Prime Plus® VET, Cartuccia calibratore con creatinina Stat Profile Prime Plus® VET, Cartucho de calibrador Stat Profile Prime Plus® VET com creatinina, Stat Profile Prime Plus® VET kalibrátor patron kreatininnel, Stat Profile Prime Plus® VET קריאטינין כיוול עם מחסנית כיוול, Stat Profile Prime Plus® VET キャリブレーター (クレアチニン), 크레아티닌 사용 Stat Profile Prime Plus® VET 보정기 카트리지, Stat Profile Prime Plus® 含肌酐的校准试剂盒

LOT 24106044

2025-10-01

| CAL A ≥ 495 mL   |              | CAL B ≥ 250 mL  |               | CAL C ≥ 245 mL   |              | CAL D ≥ 250 mL |              | CAL F ≥ 850 mL  |            | SOLN R ≥ 500 mL                 |         |
|------------------|--------------|-----------------|---------------|------------------|--------------|----------------|--------------|-----------------|------------|---------------------------------|---------|
| HCO <sub>3</sub> | 11 mmol/L    | pH              | 6.834 @ 37°C  | HCO <sub>3</sub> | 22 mmol/L    | pH             | 7.354 @ 37°C | PO <sub>2</sub> | 101.5 mmHg | KCl                             | 2 mol/L |
| Na <sup>+</sup>  | 135.5 mmol/L | H <sup>+</sup>  | 146.55 nmol/L | Cl <sup>-</sup>  | 69.0 mmol/L  | H <sup>+</sup> | 44.26 nmol/L | PO <sub>2</sub> | 13.50 kPa  | <b>WASTE CONTAINER</b> <b>W</b> |         |
| K <sup>+</sup>   | 4.0 mmol/L   | Na <sup>+</sup> | 77.0 mmol/L   | Glu              | 200 mg/dL    |                |              |                 |            |                                 |         |
| Cl <sup>-</sup>  | 108.0 mmol/L | K <sup>+</sup>  | 10.0 mmol/L   | Glu              | 11.10 mmol/L |                |              |                 |            |                                 |         |
| iCa              | 1.04 mmol/L  | iCa             | 1.92 mmol/L   | Lac              | 10.0 mmol/L  |                |              |                 |            |                                 |         |
| iCa              | 4.16 mg/dL   | iCa             | 7.68 mg/dL    | Lac              | 89.0 mg/dL   |                |              |                 |            |                                 |         |
| iMg              | 0.50 mmol/L  | iMg             | 1.50 mmol/L   | BUN              | 50 mg/dL     |                |              |                 |            |                                 |         |
| iMg              | 1.22 mg/dL   | iMg             | 3.63 mg/dL    | BUN              | 17.9 mmol/L  |                |              |                 |            |                                 |         |
| Glu              | 80 mg/dL     | tHb             | 15.4 g/dL     | Creat            | 5.0 mg/dL    |                |              |                 |            |                                 |         |
| Glu              | 4.44 mmol/L  | tHb             | 154 g/L       | Creat            | 0.44 mmol/L  |                |              |                 |            |                                 |         |
| Lac              | 2.0 mmol/L   | tHb             | 9.6 mmol/L    | Creat            | 442 μmol/L   |                |              |                 |            |                                 |         |
| Lac              | 17.8 mg/dL   |                 |               |                  |              |                |              |                 |            |                                 |         |
| BUN              | 10 mg/dL     |                 |               |                  |              |                |              |                 |            |                                 |         |
| BUN              | 3.6 mmol/L   |                 |               |                  |              |                |              |                 |            |                                 |         |
| Creat            | 1.0 mg/dL    |                 |               |                  |              |                |              |                 |            |                                 |         |
| Creat            | 0.09 mmol/L  |                 |               |                  |              |                |              |                 |            |                                 |         |
| Creat            | 88.4 μmol/L  |                 |               |                  |              |                |              |                 |            |                                 |         |

**EN**  
**Description**  
 The Cartridge contains the following aqueous solutions within individual bags: 5 Calibrators (Calibrator A, B, C, D, and F), Reference Solution, and a flexible waste bag. Each bag includes a fitter with a septum. The exposed bag fitters are arranged in a line along the rear of the cartridge box. The septa are pierced during the insertion of the cartridge into the analyzer. A charge kit consisting of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation.

**Intended Use**  
 For the quantitative determination of Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, iCa, iMg, BUN (Urea), Glucose, Lactate, Creatinine, pH, PCO<sub>2</sub>, PO<sub>2</sub>, SO<sub>2</sub> %, Hct, Hb, O<sub>2</sub>Hb, COHb, Methb, HbB, and tHb in animal whole blood using the Stat Profile Prime Plus VET Analyzer.

**Methodology**  
 Refer to Stat Profile Prime Plus VET Analyzer Instructions For Use Manual for Methodologies and Principles.

**Composition**  
 Calibrator and solution volumes are listed in the table. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF. NCCLS DOCUMENT M29-T2).

**Warnings and Cautions**  
 Refer to the Stat Profile Prime Plus VET Analyzer Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory reagents. Discard Cartridge in accordance with local guidelines. Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

**Storage**  
 Store at 15 - 30°C. Expiration date is printed on the Calibrator Cartridge.

**Directions For Use**  
**Activation of Cartridge**  
 Cal A and Cal C must be charged with creatinine syringe prior to installation of the Calibrator Cartridge. The two syringes are labeled and color coded to correspond to the fittings on the carton for Cal A and Cal C.

- Hold the syringe with tip side down and remove protective cap.
- Attach one of the enclosed needle assemblies to the syringe. Remove the protective cover from the needle.
- Match the color and label of the syringe to the appropriate fitter and insert needle.
- Slowly depress syringe plunger until the contents are dispensed.
- DO NOT PULL BACK ON THE PLUNGER TO FLUSH CONTENTS OF SYRINGE.**
- Remove needle/syringe assembly from fitter and discard in an appropriate sharps container.
- Repeat Steps 1-5 for the next Cal.
- Mix Cartridge well by gently shaking for 1 minute.

Refer to the Stat Profile Prime Plus VET Analyzer Instructions for Use Manual for complete directions. Calibrator Cartridge must be inverted gently/several times to thoroughly mix contents. Cartridge is ready for use.

**Traceability of Standards**  
 Analytes traced to NIST Standard Reference Materials wherever possible.

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