Stat Profile Prime Plus® VET Blood Gas, CO-Oximeter, Chemistry Controls Auto-Cartridge with Creatinine

Cartucho automático con creatinina para controles de química, gases en sangre y cooxímetro Stat Profile Prime Plus® VET, Cartouche automatique de contrôles de gaz du sang/CO-oxymètre, chimie Stat Profile Prime Plus® VET avec créatine, Stat Profile Prime Plus® VET Auto-Kassette mit Kreatinin für Blutgas, CO-Oximeter- und Blutchemiekontrollen, Αυτόματο φασίγγιο με κρεατινίνη επιπέδου ελέγχου χημείας αέραμας, CO-Οξυμέτρο Stat Profile Prime Plus® VET, Cartuccia con creatinina per controlli automatici chimici per gas ematici/CO-ossimetria Stat Profile Prime Plus® VET, Cartucho automático de controles de química, de CO-oxímetro e de gás no sangue Stat Profile Prime Plus® VET com creatinina, Stat Profile Prime Plus® VET vėrgaz, CO-oximėter, kėmiol kontrolė automatikais patron kreatininėl, (מחסנית בקרת כימיה אוטומטית עם קריאטינין) Stat Profile Prime Plus® VET Chemistry Controls Auto-Cartridge, Stat Profile Prime Plus® VET 血液ガス、CO オキシメーター、生化学検査用コントロール自動カートリッジ (クレアチニン)、 크라이아틴 사용 Stat Profile Prime Plus® VET 혈액 가스, CO-산소 농도계, 화학 조절액 자동 카트리지, Stat Profile Prime Plus® VET 血気、一酸化炭素血気儀、化学对照溶液自動試剤盒 (含肌酐)

LOT 24110032

CONTROL 1 2 3 4 5

2025-10-09

Table with columns for parameters (pH, H+, PCO2, PO2, SO2, Hct, Na+, K+, Cl-, iCa, Ca, iMg, Glu, Lac, BUN, Urea, Creatinine, HbF, tHb, iHb, iMg, O2Hb, COHb, MetHb, HHb, tBil, iBil) and rows for expected ranges across 5 control levels.

EN

Product Description: Material used for monitoring the performance of pH, PCO2, PO2, SO2, hematocrit (Hct), fetal hemoglobin (HbF), total hemoglobin (tHb), oxymethemoglobin (O2Hb), carboxymethemoglobin (COHb), methemoglobin (MetHb), and deoxyhemoglobin (iHb) in Levels 1, 2 and 3 as well as Na+, K+, Cl-, Ca, Mg, Glucose, Lactate, BUN (Urea), and Creatinine in Levels 4 and 5. A change kit consisting of 5 syringes filled with creatinine solution is included for cartridge activation prior to installation.

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Limitations: PCO2 values vary inversely with temperature (approximately 1%/°C). Therefore, it is critical to follow the temperature guidelines described in "Directions for Use". The Expected Range values are specific for instruments and controls manufactured by Nova Biomedical. Once installed, each Stat Profile Prime Plus VET Cartridge may be used for a maximum of 21 days from the initial installation date on the system at which time the system will indicate the cartridge is invalid.

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200 Prospect Street Waltham, MA 02454 U.S.A.