LOT 24110023

2026-04-09

CONTROL 1 2 3

15°C 30°C

IVD

Expected Ranges, Erwartungsbereiche, Αναμενόμενη περιοχή τιμών, Gamas esperadas, Plages prévues, Intervalli attesi, Intervalos esperados, Förväntade områden

Table with 4 columns: Parameter, Unit, and three Control columns (CONTROL 1, CONTROL 2, CONTROL 3) showing expected ranges for various analytes like pH, H+, PCO2, PO2, SO2, Hct, Hb, Na+, K+, Cl-, Ca++, Glu.

Product Description: An aqueous quality control material for monitoring the measurement of pH, pCO2, pO2, SO2, hematorcrit (Hct), hemoglobin (Hb), Na+, K+, Cl-, Ca+, Glucose (Glu) for use with Nova Biomedical analyzers ONLY.

Intended Use: For in vitro diagnostic use for monitoring the performance of Nova Biomedical Stat Profile pHox Plus/C Analyzers.

Methodology: Refer to Stat Profile pHox Plus/C Analyzer Instructions For Use Manual.

Composition: A buffered bicarbonate solution, each control with a known pH and known levels of Na, K, Cl, Ca+, Glu. Solutions are equilibrated with known levels of O2, CO2, and N2.

Warnings and Cautions: Must be stored at 2-28°C for at least 24 hours before opening. Intended for in vitro diagnostic use.

Storage: Store at 15-30°C; DO NOT FREEZE. Expiration date is printed on each ampule container. Do not store unused opened ampules.

Directions for use: Contents must be shaken well prior to use. Controls must be stored at approximately 24-26°C for at least 24 hours prior to opening.

Limitations: PO2 values vary inversely with temperature (approximately 1%/°C). The Expected Range values are specific for analyzers and calibrators manufactured by Nova Biomedical.

Traceability of Standards: Analytes are traced to NIST Standard Reference Materials.

Reference Intervals: Concentrations are formulated to represent three pH levels (Acidosis, Normal pH, and Alkalosis).

Users may wish to determine MEAN VALUES and EXPECTED RANGES in their own laboratory.

Expected Ranges: The EXPECTED RANGE for each analyte was determined at Nova Biomedical by using multiple runs of each level of control at 37°C on multiple analyzers.

*How to Define and Determine Reference Intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13

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