REF 48927, 48928

IVD 2°C 8°C CE

Stat Profile® pHox Ultra/Critical Care Xpress Chemistry Controls & Performance Check Solution Auto-Cartridge With Creatinine

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Table with columns for Analyte, Units, Control 4 (min-x-max), and Control 5 (min-x-max). Rows include Na+, K+, Cl-, Ca++, Mg++, Glu, Lac, BUN, Urea, Creatinine, etc.

Table with columns for PERFORMANCE CHECK and SOLN. Rows include PASS - No Action, FAIL - Replace Membrane, FEHLERFREI - Keine Maßnahme, FEHLERHAFT - Membran austauschen, etc.

Product Description: Carton contains three bags. Two bags contain assayed quality control material formulated for monitoring the measurement of Na+, K+, Cl-, Ca++, Mg+, Glucose, Lactate, BUN and Creatinine...

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Intended Use: For in vitro diagnostic use to monitor the performance of the Stat Profile pHox Ultra/Critical Care Xpress Analyzer...

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Composition: Controls (4 and 5) contain known concentrations of Na+, K+, Cl-, Ca++, Mg++, Glucose, Lactate, BUN, and Creatinine at two levels.

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Warnings and Precautions: For in vitro diagnostic use. Do not freeze. Refer to the Stat Profile pHox Ultra/Critical Care Xpress Instructions for Use Manual...

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Directions for Use: Activation of Cartridge. Each level of control must be charged with the creatinine syringe prior to installation of the Controls Auto-Cartridge.

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Limitations: This product is warranted to perform as described in this product insert, and Nova Biomedical disclaims any implied warranty of merchantability or fitness for any other purpose...

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Traceability of Controls/Solution: Analyzers are traced to NIST Standard reference materials.

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Expected Ranges: The expected range for each analyte was determined at Nova Biomedical using replicate determinations of Nova analyzers.

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References: \* NCCLS Document M29-T2. \* How to define and determine reference intervals in the clinical laboratory, approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13.

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Additional notes and footnotes regarding document references and performance metrics.

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EC REP Nova Biomedical GmbH Hessenring 13A, Geb. G 64546 Mörfelden-Walldorf, Germany