REF 52427,52861,53365,53105,53469,53470,52862,52863,53364,53466,53467,53468

# Stat Profile® Prime Calibrator Cartridge CCS/CCS Comp

Kalibratorkassette, Κασέτα Βαθμονομητή, Cartucho del calibrador, Cartouche étalon, Cartuccia calibratore, Calibrator Cartridge, Kalibratorkassett, Kalibrátorkazetta, キャリブレータカートリッジ, 측정기 카트리지, 校准液盒

LOT 24117073

2025-10-19

CAL A ≥ 500 mL		CAL B ≥ 250 mL		CAL F ≥ 720 mL		SOLN R ≥ 300 mL	
pH	7.350	pH	6.840	PCO <sub>2</sub>	51.8 mmHg	KCl	2 mol/L
PCO <sub>2</sub>	25.2 mmHg	PO <sub>2</sub>	167.0 mmHg	PO <sub>2</sub>	94.4 mmHg	WASTE CONTAINER W	
Na <sup>+</sup>	133.0 mmol/L	Na <sup>+</sup>	68.5 mmol/L				
K <sup>+</sup>	4.0 mmol/L	K <sup>+</sup>	10.0 mmol/L				
Cl <sup>-</sup>	106.0 mmol/L	Cl <sup>-</sup>	46.0 mmol/L				
iCa	1.10 mmol/L	iCa	2.20 mmol/L				
Glu	80 mg/dL	Glu	200 mg/dL				
Lac	2.0 mmol/L	Lac	10.0 mmol/L				

**EN**  
Product Description  
The above table lists the concentrations of constituents for the calibrators and solutions contained in the Calibrator Cartridge. Contains the following aqueous solutions within individual flexible bags housed in a cardboard box: 3 Calibrators (Calibrator A, B, and F), reference solution, and a flexible waste bag. Each bag includes a filament with a septa. The exposed bag filaments are arranged in a line along the rear of the cartridge box. These septa are pierced during the insertion of the cartridge into the analyzer.

Intended Use  
The Stat Profile Prime Calibrator Cartridge CCS is intended for the calibration of pH, PCO<sub>2</sub>, PO<sub>2</sub>, Hct, Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, iCa, Glucose (Glu), and Lactate (Lac) using the Stat Profile Prime CCS Analyzer.

Methodology  
Refer to Prime Analyzer Instructions For Use Manual for Test Methodologies and Principles.

Composition  
Calibrator volumes are listed in the table above. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF. NCCLS DOCUMENT M29-T2).

Warnings and Precautions  
Intended for in vitro diagnostic use.  
Refer to the Prime Analyzer Instructions for Use Manual for complete directions for use, installation instructions, recommendations for calibration intervals, troubleshooting information, Methodology and Principles of the testing procedures. Follow standard practices required for handling laboratory reagents. Discard Calibrator Cartridges and Waste products in accordance with local guide lines. Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

Storage  
Store at 15-30°C. Expiration date is printed on the front of the Calibrator Cartridge.  
Directions For Use Refer to the Prime Analyzer Instructions For Use for complete directions for installation of the cartridge, and recommendations for Calibration Intervals. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents. Cartridge is ready for use.

Traceability of Standards  
Analyses traced to NIST Standard Reference Materials.

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