Stat Profile® pHox Ultra/Critical Care Xpress Chemistry Controls

Chemiekontrollen, Χημικοί Ελεγκτές, Controles para química, Contrôles chimiques, Controlli chimici, Controles de química, Kemiska kontrollor

CONTROL 4 5

LOT 24128029

2025-05-03

LOT

CONTROL 4 24124008 2025-05-03

CONTROL 5 24124009 2025-05-03

Expected Ranges, Erwartete Bereiche, Αναμενόμενα Εύρη, Rangos esperados, Fourchettes attendues, Intervalli previsti, Gamas previstas, Förväntade områden

		CONTROL 4	CONTROL 5
Na ⁺	mmol/L	132.2 - 136.2 - 140.2	107.7 - 111.7 - 115.7
K ⁺	mmol/L	3.66 - 3.91 - 4.16	5.95 - 6.25 - 6.55
Cl ⁻	mmol/L	119.1 - 123.6 - 128.1	93.0 - 97.5 - 102.0
Ca ⁺⁺	mmol/L	0.90 - 0.98 - 1.06	1.38 - 1.48 - 1.62
Mg ⁺⁺	mmol/L	0.53 - 0.60 - 0.67	0.95 - 1.10 - 1.25
Glu	mg/dL	71 - 78 - 85	251 - 276 - 301
Glu	mmol/L	3.94 - 4.33 - 4.72	13.93 - 15.32 - 16.71
Lac	mmol/L	1.7 - 2.0 - 2.3	6.0 - 6.7 - 7.4
BUN	mg/dL	14 - 19 - 24	51 - 58 - 65
BUN	mmol/L	5.0 - 6.8 - 8.6	18.2 - 20.7 - 23.2
Urea	mg/dL	30.0 - 40.8 - 51.5	109.4 - 124.4 - 139.4
Urea	mmol/L	5.0 - 6.8 - 8.6	18.2 - 20.7 - 23.2
Creatinine	mg/dL	0.8 - 1.1 - 1.4	5.5 - 6.4 - 7.3
Creatinine	μmol/L	70.8 - 97.3 - 123.9	486.6 - 566.2 - 645.8

Product Description
Contains aqueous quality control material formulated at two levels: Control 4 and 5, for monitoring the performance of Na⁺, K⁺, Cl⁻, Ca⁺⁺, Mg⁺⁺, Glucose, Lactate, BUN (urea), and Creatinine on the Stat Profile pHox Ultra/Critical Care Xpress Analyzer.

Intended Use
For in vitro diagnostic use to monitor the performance of the Stat Profile pHox Ultra/Critical Care Xpress Analyzer.

Methodology
Refer to Stat Profile pHox Ultra/Critical Care Xpress Analyzer Instruction for Use Manual for Methodology and Principles.

Composition
Chemistry Control 4 and 5 contain known concentrations of Na⁺, K⁺, Cl⁻, Ca⁺⁺, Mg⁺⁺, Glucose, Lactate, BUN, and Creatinine in a buffered solution. Each ampule contains a minimum volume of 1.7 mL. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials¹.

Warning and Precautions
For in vitro diagnostic use. Store at 2-8°C. Do not freeze. Refer to the Stat Profile pHox Ultra/Critical Care Xpress instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory reagents.

Storage
Store at 2-8°C. Do not freeze. Each ampule has a lot number and expiration date printed on the label.

Directions for Use
Shake ampule before opening. Snap open ampule (protecting fingers with gauze or gloves), and within ONE MINUTE, aspirate solution with sample probe. See Instructions for Use Manual for complete directions.

Limitations
For use on Nova Biomedical Analyzers Only. Performance characteristics have not been established for use on other manufacturer's analyzers.

This product is warranted to perform as described in this product insert, and Nova Biomedical disclaims any implied warranty of merchantability or fitness for any other purpose, and in no event shall Nova Biomedical be liable for any consequential damages arising out of the aforesaid expressed warranty.

Traceability of Standards
Analytes are traced to NIST Standard reference materials.

Reference Intervals
Concentrations are formulated at normal and abnormal expected values in patient blood. The expected clinical range of these values in patient blood is referenced in Tietz, NW ed. 1985 Textbook of Clinical Chemistry, W.B. Saunders Co. Users may wish to determine Mean values and expected ranges in their own laboratory.²

Expected ranges
The expected range for each analyte was determined at Nova Biomedical using replicate determinations on Nova analyzers. The expected range indicates the maximum deviations from the mean value that may be expected under differing laboratory conditions for instruments operating within specifications. Refer to Expected Ranges Table.

¹ NCCLS Document M29-T2.
² How to define and determine reference intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13.

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