REF 52714

Stat Profile® Prime Ampuled Control ABG/CCS CONTROL 1 2 3

LOT

24123023 · 2025-10-30

Rx Only

IVD

15°C - 30°C

EN

Product Description
An aqueous quality control material for monitoring the measurement of pH, PCO₂, PO₂, Hct, Na, K, Cl, iCa, Glucose (Glu), and Lactate (Lac) for use with Nova Biomedical analyzers ONLY. Formulated at 3 levels:

CONTROL 1 Acidosis, with High Electrolyte, Low Normal Glucose, Normal Lactate
CONTROL 2 Normal pH, with Normal Electrolyte, High Glucose, High Lactate
CONTROL 3 Alkalosis, with Low Electrolyte, High Abnormal Glucose, High Abnormal Lactate

Intended Use
The StatProfile Prime Ampuled Control ABG/CCS is a Quality Control material intended for *in vitro* diagnostic use by healthcare professionals for monitoring the performance of the StatProfile Prime CCS Analyzer.
Methodology
Refer to Prime Analyzer Instructions For Use Manual.

Composition
A buffered bicarbonate solution, each control with a known pH and known levels of Na, K, Cl, iCa, Glu, and Lac. Solutions are equilibrated with known levels of O₂, CO₂, and N₂. The reflectance characteristics give a signal that is equivalent to a known hematocrit value in whole blood. Mold inhibited. Each ampule contains 1.7 ml volume. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF. NCCLS DOCUMENT M29-T2)

Warnings and Cautions:
Must be stored at 24 - 26°C for at least 24 hours before opening. Intended for *in vitro* diagnostic use. Refer to the Prime Analyzer Instructions for Use for complete directions for use, including application of the control sample onto the analyzer, recommendations for use of controls, troubleshooting information and methodology and principles of the testing procedures. Follow standard practices required for handling laboratory reagents. Once ampule is opened, discard unused portion in accordance with local guide lines.

Storage
Store at 15-30°C; DO NOT FREEZE. Expiration date is printed on each ampule container. Do not store unused opened ampules.

Directions for use
Contents must be shaken well prior to use. Controls must be stored at approximately 24-26°C for at least 24 hours prior to opening. Refer to Analyzer Instructions for Use for complete instructions. Verify that the lot number appearing on the Expected Ranges Table is the same lot number indicated on the control ampule.

Limitations
PO₂ values vary inversely with temperature (approximately 1%/°C).
The Expected Range values are specific for instruments and controls manufactured by Nova Biomedical.

Traceability of Standards
Analytes are traced to NIST Standard Reference Materials.

Reference Intervals
Concentrations are formulated to represent 3 pH levels (Acidosis, Normal pH, and Alkalosis). The expected clinical range of values for these analytes in patient blood is referenced in Tietz, N.W. ed 1986 Textbook of Clinical Chemistry, W.B. Saunders Co. Users may wish to determine MEAN VALUES and EXPECTED RANGES in their own laboratory.

Expected Ranges The EXPECTED RANGE for each analyte was determined at Nova Biomedical by using multiple runs of each level of control at 37°C on multiple instruments. The EXPECTED RANGE indicates the maximum deviations from the mean value that may be expected under differing laboratory conditions for instruments operating within specifications.

Refer to Expected Ranges Table
Conversion Factors: For Ca⁺⁺, Glucose, and Lactate
Ca⁺⁺ in mmol/L x 4 = mg/dL
Glucose in mg/dL ÷ 18.016 = mmol/L
Lactate in mmol/L x 8.9 = mg/dL

How to Define and Determine Reference Intervals in the clinical laboratory, approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13.

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EC REP Nova Biomedical GmbH Hessenring 13A, Geb. G 64546 Mörfelden-Walldorf Germany
Waltham, MA 02454 U.S.A.