REF 48905, 48906

IVD 2°C - 8°C CE

# Stat Profile® pHox Ultra/Critical Care Xpress Chemistry Controls & Performance Check Solution Auto-Cartridge

Chemiekontrollen und Leistungstest-Lösung Auto-QK-Kassette, Αυτόματη Κασέτα Χημικών Ελεγκτών & Διαλύματος Ελέγχου Απόδοσης, Cartucho automático para Controles para química y Solución de verificación de comportamiento, Cartouche CQ de contrôles chimiques et de solution de contrôle de linéarité, Cartuccia automatica controlli chimici e soluzione controllo prestazioni, Controlos de química & Cartucho automático de Solução de Verificação de Desempenho, Autokassett med kemiska kontroll och prestandakontrollösning

LOT 24138012 2025-11-07

CONTROL 4 5 PERFORMANCE CHECK SOLN

### Expected Ranges, Erwartete Bereiche, Αναμενόμενα Εύρη, Rangos esperados, Fourchettes attendues, Intervalli previsti, Gamas previstas, Förväntade områden

		CONTROL 4	CONTROL 5
		min - $\bar{x}$ - max	min - $\bar{x}$ - max
Na <sup>+</sup>	mmol/L	139.2 - 143.2 - 147.2	112.8 - 116.8 - 120.8
K <sup>+</sup>	mmol/L	3.80 - 4.05 - 4.30	6.28 - 6.58 - 6.88
Cl <sup>-</sup>	mmol/L	122.2 - 126.7 - 131.2	91.5 - 96.0 - 100.5
Ca <sup>++</sup>	mmol/L	1.03 - 1.11 - 1.19	1.52 - 1.62 - 1.76
Mg <sup>++</sup>	mmol/L	0.54 - 0.63 - 0.72	0.98 - 1.13 - 1.28
Glu	mg/dL	73 - 80 - 87	255 - 280 - 305
Glu	mmol/L	4.05 - 4.44 - 4.83	14.15 - 15.54 - 16.93
Lac	mmol/L	1.7 - 2.0 - 2.3	6.3 - 7.0 - 7.7
BUN	mg/dL	14 - 19 - 24	49 - 56 - 63
BUN	mmol/L	5.0 - 6.8 - 8.6	17.5 - 20.0 - 22.5
Urea	mg/dl	30.0 - 40.8 - 51.5	105.1 - 120.1 - 135.1
Urea	mmol/L	5.0 - 6.8 - 8.6	17.5 - 20.0 - 22.5

PERFORMANCE CHECK	SOLN
PASS - No Action	
FAIL - Replace Membrane	
FEHLERFREI - Keine Maßnahme	
FEHLERHAFT - Membran austauschen	
ΕΠΙΤΥΧΙΑ - Καμία Ενέργεια	
ΑΠΟΤΥΧΙΑ - Αντικαταστήστε τη μεμβράνη	
GLUCOSE	
PASO - Ninguna acción	
NO PASO - Sustituir membrana	
SUCES - Aucune action	
ECHEC - Remplacer la membrane	
RIUSCITO - Nessuna azione	
NON RIUSCITO - Sostituisci membrana	
PASSO - Nenhuma acción	
FALHOU - Substituir membrana	
OK - Ingen åtgärd	
FEL - Byt ut membran	

**Product Description**  
Carton contains three bags. Two bags contain assayed quality control material formulated for monitoring the measurement of Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, Ca<sup>++</sup>, Mg<sup>++</sup>, Glucose, Lactate, and BUN (Controls 4 and 5). One bag contains an aqueous quality control material to verify glucose membrane performance (Performance Check Solution, Control 6).

**Intended Use**  
For in vitro diagnostic use to monitor the performance of the Stat Profile pHox Ultra/Critical Care Xpress Analyzer and to insure optimum performance of glucose membrane. Performance Check solution should be analyzed before controls are run and after changing glucose membrane.

**Methodology**  
Refer to Stat Profile pHox Ultra/Critical Care Xpress Analyzer Instruction for Use Manual.

**Composition**  
Controls (Controls 4 and 5) contain known concentrations of Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, Ca<sup>++</sup>, Mg<sup>++</sup>, Glucose, Lactate, and BUN at two levels. Performance Check Solution (Control 6) is a buffered salt solution containing Glucose (1000mg/dL). Each bag contains a minimum volume of 90 mL. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials.

**Warnings and Precautions**  
For in vitro diagnostic use. Do not freeze. Refer to the Stat Profile pHox Ultra/Critical Care Xpress Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory controls.

**Storage**  
Store at 2-8°C. Do not freeze. Each cartridge has lot number and expiration date printed on the label.

**Directions for Use**  
Refer to Analyzer Instructions for Use Manual for complete directions. Performance Check solution should be analyzed before controls are run and after changing glucose membranes. Verify that the lot number appearing on the expected ranges table is the same as the lot number printed on the cartridge.

**Limitations**  
For use on Nova Biomedical Analyzers Only. Performance characteristics have not been established for use on other manufacturer's analyzers. Once installed, each Control Auto-Cartridge may be used for a maximum of 35 days from the date of initial activation and installation. The Analyzer will indicate that the cartridge is invalid at that time. This product is warranted to perform as described in this product insert, and Nova Biomedical disclaims any implied warranty of merchantability or fitness for any other purpose, and in no event shall Nova Biomedical be liable for any consequential damages arising out of the aforesaid expressed warranty.

**Traceability of Standards**  
Analyses are traced to NIST Standard reference materials.

**Reference Intervals**  
Concentrations are formulated at normal and abnormal levels expected in patient blood. The expected clinical ranges of patient blood is referenced in Tietz, N.W. ed. 1996 Textbook of Clinical Chemistry, W.B. Saunders Co. Users may wish to determine Expected Ranges in their own laboratory.<sup>1</sup>

**Expected Ranges**  
The expected range for each analyte was determined at Nova Biomedical using replicate determinations on Nova analyzers. The expected range indicates the maximum deviations from the mean value that may be expected under differing laboratory conditions for analyzers operating within specifications. Refer to Expected Ranges table for Control 4 and 5 values. Performance Check Solution results for glucose will be Pass or Fail.

<sup>1</sup> NCCLS Document M29-T2.  
<sup>2</sup> How to define and determine reference intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13.

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## Tracing of the standards

The analyses are traced according to the parameters of the reference standard of the NIST (National Institute of Standards and Technology).

## Reference intervals

The concentrations are formulated at the normal and abnormal expected levels in the blood of patients.

The expected clinical values measured in the blood of patients are mentioned in Tietz, N.W. ed. 1986. Textbook of Clinical Chemistry, W.B. Saunders Co.

Users may want to determine the expected intervals in their own laboratories.<sup>2</sup>

## Expected ranges

The expected range for each analysis was determined by Nova Biomedical, repeating several times each level on Nova analyzers. The expected range indicates the maximum deviations from the mean value expected in different laboratory controls for analyzers that function within the specifications. Consult the table of Expected Ranges for the values of control 4 and 5. The results of the linearity solution for glucose will be Success (Pass) or Failure (Fail).

## Traceability of standards

The analyses are traced according to the parameters of the reference standard of the NIST (National Institute of Standards and Technology).

## Reference intervals

The concentrations are formulated at the normal and abnormal expected levels in the blood of patients.

The expected clinical values measured in the blood of patients are mentioned in Tietz, N.W. ed. 1986. Textbook of Clinical Chemistry, W.B. Saunders Co.

Users may want to determine the expected intervals in their own laboratories.<sup>2</sup>

## Expected ranges

The expected range for each analysis was determined by Nova Biomedical, repeating several times each level on Nova analyzers. The expected range indicates the maximum deviations from the mean value expected in different laboratory controls for analyzers that function within the specifications. Consult the table of Expected Ranges for the values of control 4 and 5. The results of the linearity solution for glucose will be Success (Pass) or Failure (Fail).

## Traceability of standards

The analyses are traced according to the parameters of the reference standard of the NIST (National Institute of Standards and Technology).

## Reference intervals

The concentrations are formulated at the normal and abnormal expected levels in the blood of patients.

The expected clinical values measured in the blood of patients are mentioned in Tietz, N.W. ed. 1986. Textbook of Clinical Chemistry, W.B. Saunders Co.

Users may want to determine the expected intervals in their own laboratories.<sup>2</sup>

## Expected ranges

The expected range for each analysis was determined by Nova Biomedical, repeating several times each level on Nova analyzers. The expected range indicates the maximum deviations from the mean value expected in different laboratory controls for analyzers that function within the specifications. Consult the table of Expected Ranges for the values of control 4 and 5. The results of the linearity solution for glucose will be Success (Pass) or Failure (Fail).

## Traceability of standards

The analyses are traced according to the parameters of the reference standard of the NIST (National Institute of Standards and Technology).

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