Stat Profile Prime Plus® Blood Gas, CO-Oximeter, Chemistry Controls Auto-Cartridge with Creatinine

Cartucho automático con creatinina para controles de química, gases en sangre y cooxímetro Stat Profile Prime Plus®. Cartouche automatique de contrôles de gaz du sang/CO-oxymètre, chimie

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Table with 6 columns: Parameter, CONTROL 1, CONTROL 2, CONTROL 3, CONTROL 4, CONTROL 5. Rows include pH, H+, PCO2, PO2, SO2, Hct, Na+, K+, Cl-, iCa, iMg, Glu, Lac, Bun, Urea, Creatinine, HbF, tHb, tHb, tHb, O2Hb, COHb, MetHb, HHb, tBil, tBil, tBil.

EN

Product Description: Analytical quality control material for monitoring the performance of pH, PCO2, PO2, SO2, hematocrit (Hct), fetal hemoglobin (HbF), total hemoglobin (tHb), total bilirubin (tBil), oxyhemoglobin (O2Hb), carboxyhemoglobin (COHb), methemoglobin (MetHb), and desoxyhemoglobin (dHb) in Levels 1, 2 and 3 as well as Na+, K+, Cl-, Ca, Mg, glucose, lactate, BUN (urea), and Creatinine in Level 4 and 5. A charge kit consisting of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation.

Directions for use: Ensure controls are room temperature prior to installation. Level 4 and Level 5 must be charged with creatinine syringes prior to installation of the Control Cartridge. The 2 syringes are labeled and color coded to correspond to the terminals on the cartridge. Activate the cartridge as follows: 1. Hold the syringe with top side down and remove protective cap. 2. Attach one of the enclosed needle assemblies to the syringe. Remove the protective cover from the needle. 3. Moisten the color and label of the syringe to the appropriate ferment and insert needle. 4. Slowly depress syringe plunger until the contents are dispensed. DO NOT PULL BACK ON THE PLUNGER TO FLUSH THE CONTENTS OF THE SYRINGE. 5. Remove needle/syringe assembly from ferment and discard in an appropriate sharps container. 6. Repeat Steps 1-4 for the next control. 7. Mix Cartridge well by gently swirling for 1 minute. Cartridge is ready for use.

Limitations: PO2 values vary inversely with temperature (approximately 1% / °C). Therefore, it is critical to follow the temperature guidelines described in "Directions for Use". The Expected Range values are specific for instruments and controls manufactured by Nova Biomedical. Contact installed, each Stat Profile Prime Plus Cartridge may be used for a maximum of 21 days from the initial installation date on the system at which time the system will indicate the cartridge is invalid. Each cartridge may be inserted and removed from the analyzer a maximum of 5 times.

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