Stat Profile Prime Plus® VET Blood Gas, CO-Oximeter, Chemistry Controls Auto-Cartridge with Creatinine

Cartucho automático con creatinina para controles de química, gases en sangre y cooxímetro Stat Profile Prime Plus® VET, Cartouche automatique de contrôles de gaz du sang/CO-oxymètre, chimio Stat Profile Prime Plus® VET avec créatine, Stat Profile Prime Plus® VET Auto-Kassette mit Kreatinin für Blutgas-, CO-Oximeter- und Blutchemiekontrollen, Αυτόματη συσκευή με κρεατινίνη επιπέδου ελάγχου...

LOT 24158023

CONTROL 1 2 3 4 5

2025-11-23

Table with 7 columns: Parameter, CONTROL 1, CONTROL 2, CONTROL 3, CONTROL 4, CONTROL 5. Rows include pH, H+, PCO2, PO2, SO2, Hct, Na+, K+, Cl-, iCa, iCa, iMg, iMg, Glu, Glu, Lac, Lac, BUN, BUN, Urea, Urea, Creatinine, Creatinine, Creatinine, HbF, iHb, iHb, iHb, O2Hb, COHb, MetHb, HHb, iBil, iBil, iBil.

EN

Product Description: Aqueous quality control material for monitoring the performance of pH, PCO2, PO2, SO2, hematocrit (Hct), fetal hemoglobin (HbF), total hemoglobin (HbT), oxymethemoglobin (COHb), methemoglobin (MetHb), and deoxyhemoglobin (HbD) in Levels 1, 2 and 3 as well as Na+, K+, Cl-, Ca, Mg, Glucose, Lactate, BUN (Urea), and Creatinine in Levels 4 and 5. A cartridge kit consisting of 5 syringes filled with creatinine solution is included for cartridge activation prior to installation.

Warnings and Cautions: DO NOT FREEZE. Mix the cartridge by gently inverting for several seconds. DO NOT SHAKE CARTRIDGE. Refer to Stat Profile Prime Plus VET Analyzer Instructions for Use Manual for complete information.

Storage: Store at 2-8°C (37-46°F). DO NOT FREEZE.

Directions for use: Ensure controls are room temperature prior to installation. Level 4 and Level 5 must be charged with creatinine syringes prior to installation of the Calibrator Cartridge. The 2 syringes are labeled and color coded to correspond to the openings on the cartridge.

Traceability of Standards: Total hemoglobin (HbT) and methemoglobin (MetHb) are traceable by using Cyanmethemoglobin method. Carboxyhemoglobin (COHb) and Oxymethemoglobin (O2Hb) are traceable using Spectrophotometry. Analytes are traced to NIST Standard Reference Materials.

Expected Ranges: The expected range for each parameter was determined at Nova Biomedical using replicate determinations on Nova analyzers. The expected range indicates the maximum deviations from the Mean Value that may be expected after allowing laboratory conditions for instrument operation within specifications.

INCCLS Document M29-12: Show to Define and Determine Reference Intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13.

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