REF 58395,58396

# Stat Profile Prime Plus® VET Calibrator Cartridge

Stat Profile Prime Plus® VET Καλιβριετή, Φυσιόγιο βαθμονόμησης VET Stat Profile Prime Plus®, Cartucho de calibrador Stat Profile Prime Plus® VET, Cartouche de calibreurs Stat Profile Prime Plus® VET, Cartuccia calibratore VET Stat Profile Prime Plus®, Cartucho de calibrador VET Stat Profile Prime Plus®, Stat Profile Prime Plus® VET kalibrátor patron, Stat Profile Prime Plus® VET מכונת כיול, Stat Profile Prime Plus® VET キャリブレーターカートリッジ, Stat Profile Prime Plus® VET 보정기 카트리지, Stat Profile Prime Plus® VET 校正试剂盒

LOT 24180039

2025-12-19

| CAL A            | ≥ 495 mL     | CAL B           | ≥ 250 mL      | CAL C            | ≥ 245 mL     | CAL D          | ≥ 250 mL     | CAL F           | ≥ 850 mL   | SOLN R | ≥ 500 mL |
|------------------|--------------|-----------------|---------------|------------------|--------------|----------------|--------------|-----------------|------------|--------|----------|
| HCO <sub>3</sub> | 11 mmol/L    | pH              | 6.834 @ 37°C  | HCO <sub>3</sub> | 22 mmol/L    | pH             | 7.354 @ 37°C | PO <sub>2</sub> | 103.0 mmHg | KCl    | 2 mol/L  |
| Na <sup>+</sup>  | 135.5 mmol/L | H <sup>+</sup>  | 146.55 nmol/L | Cl <sup>-</sup>  | 69.0 mmol/L  | H <sup>+</sup> | 44.26 nmol/L | PO <sub>2</sub> | 13.70 kPa  |        |          |
| K <sup>+</sup>   | 4.0 mmol/L   | Na <sup>+</sup> | 77.0 mmol/L   | Glu              | 200 mg/dL    |                |              |                 |            |        |          |
| Cl <sup>-</sup>  | 108.0 mmol/L | K <sup>+</sup>  | 10.0 mmol/L   | Glu              | 11.10 mmol/L |                |              |                 |            |        |          |
| iCa              | 1.04 mmol/L  | iCa             | 1.92 mmol/L   | Lac              | 10.0 mmol/L  |                |              |                 |            |        |          |
| iCa              | 4.16 mg/dL   | iCa             | 7.68 mg/dL    | Lac              | 89.0 mg/dL   |                |              |                 |            |        |          |
| iMg              | 0.50 mmol/L  | iMg             | 1.50 mmol/L   | BUN              | 50 mg/dL     |                |              |                 |            |        |          |
| iMg              | 1.22 mg/dL   | iMg             | 3.63 mg/dL    | BUN              | 17.9 mmol/L  |                |              |                 |            |        |          |
| Glu              | 80 mg/dL     | tHb             | 15.4 g/dL     |                  |              |                |              |                 |            |        |          |
| Glu              | 4.44 mmol/L  | tHb             | 154 g/L       |                  |              |                |              |                 |            |        |          |
| Lac              | 2.0 mmol/L   | tHb             | 9.6 mmol/L    |                  |              |                |              |                 |            |        |          |
| Lac              | 17.8 mg/dL   |                 |               |                  |              |                |              |                 |            |        |          |
| BUN              | 10 mg/dL     |                 |               |                  |              |                |              |                 |            |        |          |
| BUN              | 3.6 mmol/L   |                 |               |                  |              |                |              |                 |            |        |          |

**EN**  
**Description**  
 The Cartridge contains the following aqueous solutions within individual bags: 5 Calibrators (Calibrator A, B, C, D, and F), Reference Solution, and a flexible waste bag. Each bag includes a filament with a septum. The exposed bag filaments are arranged in a line along the rear of the cartridge box. The septa are pierced during the insertion of the cartridge into the analyzer. A charge kit consisting of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation.  
**Intended Use**  
 For the quantitative determination of Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, iCa, iMg, BUN (Urea), Glucose, Lactate, pH, PCO<sub>2</sub>, PO<sub>2</sub>, SO<sub>2</sub>, Hct, tHb, O<sub>2</sub>Hb, COHb, MetHb, Hb, and tBil in animal whole blood using the Stat Profile Prime Plus VET Analyzer.  
**Methodology**  
 Refer to Stat Profile Prime Plus VET Analyzer Instructions For Use Manual for Methodologies and Principles.  
**Composition**  
 Calibrator and solution volumes are listed in the table. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF. NCCLS DOCUMENT M29-T2).  
**Warnings and Cautions**  
 Refer to the Stat Profile Prime Plus VET Analyzer Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory reagents. Discard Cartridge in accordance with local guidelines. Waste Bag Contains Preservative (eye, skin, and respiratory irritant).  
**Storage**  
 Store at 15 - 30°C. Expiration date is printed on the Calibrator Cartridge.  
**Directions For Use**  
 Refer to the Stat Profile Prime Plus VET Analyzer Instructions for Use for complete directions. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents.  
**Traceability of Standards**  
 Analytes traced to NIST Standard Reference Materials wherever possible.

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