
REF 48927, 48928

Stat Profile® pH0x Ultra/Critical Care Xpress Chemistry Controls & Performance Check Solution Auto-Cartridge With Creatinine

Chemiekontrollen und Leistungstest-Lösung Auto-QK-Kassette mit Creatinin, Αυτόματη Κασέτα Χημικών Ελεγκτών & Διαλύματος Ελέγχου Απόδοσης Με Κρεατινίνη, Cartucho automático para Controles para química y Solución de verificación de comportamiento con creatinina, Cartouche CQ de contrôles chimiques et de solution de contrôle de linéarité avec créatinine, Cartuccia automatica controlli chimici e soluzione controllo prestazioni con creatinina, Controlos de química & Cartucho automático de Solução de Verificação de Desempenho com creatinina, Autokassett med kemiska kontrollor och prestandakontrollösning med kreatinin

LOT 24169078 2025-12-05 CONTROL 4 5 PERFORMANCE CHECK SOLN

Expected Ranges, Erwartete Bereiche, Αναμενόμενα Εύρη, Rangos esperados, Fourchettes attendues, Intervalli previsti, Gamas previstas, Förväntade områden

		CONTROL 4		CONTROL 5			
		min	\bar{x}	max	min	\bar{x}	max
Na ⁺	mmol/L	136.0	140.0	144.0	110.5	114.5	118.5
K ⁺	mmol/L	3.72	3.97	4.22	6.10	6.40	6.70
Cl ⁻	mmol/L	121.3	125.8	130.3	91.5	96.0	100.5
Ca ⁺⁺	mmol/L	1.03	1.11	1.19	1.51	1.63	1.75
Mg ⁺⁺	mmol/L	0.55	0.64	0.73	1.00	1.15	1.30
Glu	mg/dL	73	80	87	255	280	305
Glu	mmol/L	4.05	4.44	4.83	14.15	15.54	16.93
Lac	mmol/L	1.7	2.0	2.3	6.3	7.0	7.7
BUN	mg/dL	14	19	24	49	56	63
BUN	mmol/L	5.0	6.8	8.6	17.5	20.0	22.5
Urea	mg/dL	30.0	40.8	51.5	105.1	120.1	135.1
Urea	mmol/L	5.0	6.8	8.6	17.5	20.0	22.5
Creatinine	mg/dL	0.7	1.0	1.3	5.2	6.2	7.2
Creatinine	μmol/L	61.9	88.5	115.0	460.0	548.5	637.0

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Product Description: Carton contains three bags. Two bags contain assayed quality control material formulated for monitoring the measurement of Na⁺, K⁺, Cl⁻, Ca⁺⁺, Mg⁺⁺, Glucose, Lactate, BUN and Creatinine (Control 4 and 5). One bag contains an aqueous quality control material to verify glucose and creatinine membrane performance (Performance Check Solution).

Intended Use: For in vitro diagnostic use to monitor the performance of the Stat Profile pH0x Ultra/Critical Care Xpress Analyzer and to ensure optimum performance of glucose and creatinine membranes.

Methodology: Refer to Stat Profile pH0x Ultra/Critical Care Xpress Analyzer Instructions for Use Manual.

Composition: Controls (4 and 5) contain known concentrations of Na⁺, K⁺, Cl⁻, Ca⁺⁺, Mg⁺⁺, Glucose, Lactate, BUN, and Creatinine at two levels. Performance Check Solution is a buffered salt solution containing known levels of Glucose, Creatinine, and Creatine.

Warnings and Precautions: For in vitro diagnostic use. Do not freeze. Refer to the Stat Profile pH0x Ultra/Critical Care Xpress Instructions for Use Manual for complete directions.

Storage: Store at 2-8°C. Do not freeze. Each cartridge has lot number and expiration date printed on the label.

Directions for Use: Activation of Cartridge. Each level of control must be charged with the creatinine syringe prior to installation of the Controls Auto-Cartridge. The three syringes are labeled and color coded to correspond to the items on the carton for each level of control.

Limitations: For use on Nova Biomedical Analyzers Only. Once installed, the Auto-Cartridge may be used for a maximum of 21 days from the date of initial activation and installation.

This product is warranted to perform as described in this product insert, and Nova Biomedical declines any implied warranty of merchantability or fitness for any other purpose.

Traceability of Controls/Solution: Analyzers are traced to NIST Standard reference materials.

Reference Intervals: Concentrations are formulated at normal and abnormal levels expected in patient blood. The expected clinical ranges of patient blood is referenced in Tietz, N.W. ed. 1986 Textbook of Clinical Chemistry, W.B. Saunders Co.

Expected Ranges: The expected range for each analyte was determined at Nova Biomedical using replicate determinations on Nova analyzers. The expected range indicates the maximum deviations from the mean value that may be expected under differing laboratory conditions for analyzers operating within specifications.

NCLS Document M29-72: How to define and determine reference intervals in the clinical laboratory; approved guideline-second edition, NCLS C28-A2, Volume 20, Number 13.

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