REF 57831, 57832, 57833, 57834, 57835, 65142

Rx Only IVD 15°C - 30°C CE

Stat Profile Prime Plus® Calibrator Cartridge With Creatinine

Stat Profile Prime Plus® Καλιβριεστική με Κρεατίνη, Stat Profile Prime Plus® Φυσιόγιο μέσω βαθμονόμησης με κρεατίνη, Cartucho de calibrador con creatinina Stat Profile Prime Plus®, Cartouche de calibrateurs Stat Profile Prime Plus® avec créatinine, Cartuccia calibratore con creatinina Stat Profile Prime Plus®, Cartucho de calibrador Stat Profile Prime Plus® com creatinina, Stat Profile Prime Plus® kalibrátor patron kreatininne, סט פרופיל פרים פלוס עם קריאטיν, Stat Profile Prime Plus®, Stat Profile Prime Plus® キャリブレーションカートリッジ (クレアチニン), 크레아티닌 사용 Stat Profile Prime Plus® 보정기 카트리지, Stat Profile Prime Plus® 含肌酐的校准试剂盒

LOT 24190024

2025-11-28

| CAL A | | CAL B | | CAL C | | CAL D | | CAL F | | SOLN R | |
|------------------|--------------|-----------------|---------------|------------------|--------------|----------------|--------------|-----------------|------------|---|---------|
| HCO ₃ | 11 mmol/L | pH | 6.834 @ 37°C | HCO ₃ | 22 mmol/L | pH | 7.354 @ 37°C | PO ₂ | 103.5 mmHg | KCl | 2 mol/L |
| Na ⁺ | 135.5 mmol/L | H ⁺ | 146.55 nmol/L | Cl ⁻ | 69.0 mmol/L | H ⁺ | 44.26 nmol/L | PO ₂ | 13.77 kPa |  | |
| K ⁺ | 4.0 mmol/L | Na ⁺ | 77.0 mmol/L | Glu | 200 mg/dL | | | | | | |
| Cl ⁻ | 108.0 mmol/L | K ⁺ | 10.0 mmol/L | Glu | 11.10 mmol/L | | | | | | |
| iCa | 1.04 mmol/L | iCa | 1.92 mmol/L | Lac | 10.0 mmol/L | | | | | | |
| iCa | 4.16 mg/dL | iCa | 7.68 mg/dL | Lac | 89.0 mg/dL | | | | | | |
| iMg | 0.50 mmol/L | iMg | 1.50 mmol/L | BUN | 50 mg/dL | | | | | | |
| iMg | 1.22 mg/dL | iMg | 3.63 mg/dL | BUN | 17.9 mmol/L | | | | | | |
| Glu | 80 mg/dL | tHb | 15.4 g/dL | Creat | 5.0 mg/dL | | | | | | |
| Glu | 4.44 mmol/L | tHb | 154 g/L | Creat | 0.44 mmol/L | | | | | | |
| Lac | 2.0 mmol/L | tHb | 9.6 mmol/L | Creat | 442 μmol/L | | | | | | |
| Lac | 17.8 mg/dL | | | | | | | | | | |
| BUN | 10 mg/dL | | | | | | | | | | |
| BUN | 3.6 mmol/L | | | | | | | | | | |
| Creat | 1.0 mg/dL | | | | | | | | | | |
| Creat | 0.09 mmol/L | | | | | | | | | | |
| Creat | 88.4 μmol/L | | | | | | | | | | |

EN
Description
 The Cartridge contains the following aqueous solutions within individual bags: 5 Calibrators (Calibrator A, B, C, D, and F), Reference Solution, and a flexible waste bag. Each bag includes a filament with a septum. The exposed bag filaments are arranged in a line along the rear of the cartridge box. The septa are pierced during the insertion of the cartridge into the analyzer. A charge kit consisting of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation.

Intended Use
 For the quantitative determination of Na⁺, K⁺, Cl⁻, iCa, iMg, BUN (Urea), Glucose, Lactate, Creatinine, pH, PCO₂, PO₂, SO₂, Hct, Hb, HbF, O₂Hb, COHb, MetHb, HbB, and IBI[†] in human blood using the Stat Profile Prime Plus Analyzer.

Methodology
 Refer to Stat Profile Prime Plus Analyzer Instructions For Use Manual for Methodologies and Principles.

Composition
 The calibrators are aqueous salt solutions containing analytes at concentrations listed in the table above. Calibrator B also contains dye. Contain no constituents of human origin, however, good laboratory practice should be followed during the handling of these materials. (REF. NCCLS DOCUMENT M29-T2).

Warnings and Cautions
 Intended for *in vitro* diagnostic use.
 Refer to the Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory reagents. Discard Cartridge in accordance with local guidelines.
 Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

Storage
 Store at 15 - 30°C. Expiration date is printed on the Calibrator Cartridge.

Directions For Use
Activation of Cartridge
 Cal A and Cal C must be charged with creatinine syringe prior to installation of the Calibrator Cartridge. The two syringes are labeled and color coded to correspond to the filaments on the carton for Cal A and Cal C.
 1. Hold the syringe with tip side down and remove protective cap.
 2. Attach one of the enclosed needle assemblies to the syringe. Remove the protective cover from the needle.
 3. Match the color and label of the syringe to the appropriate filament and insert needle.
 4. Slowly depress syringe plunger until the contents are dispensed.
DO NOT PULL BACK ON THE PLUNGER TO FLUSH CONTENTS OF SYRINGE.
 5. Remove needle/syringe assembly from filament and discard in an appropriate sharps container.
 6. Repeat Steps 1-5 for the next Cal.
 7. Mix Cartridge well by gently shaking for 1 minute.

Refer to the Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete directions. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents. Cartridge is ready for use.

*Not available in the USA or for Point-of-Care/Near-Patient Testing use.

Traceability of Standards
 Analytes traced to NIST Standard Reference Materials wherever possible.

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EC REP Nova Biomedical GmbH Hessenring 13A, Geb. G 64546 Mörfelden-Walldorf Germany
Waltham, MA 02454 U.S.A.