Stat Profile Prime Plus® VET Blood Gas, CO-Oximeter, Chemistry Controls Auto-Cartridge with Creatinine

Cartucho automático con creatinina para controles de química, gases en sangre y cooxímetro Stat Profile Prime Plus® VET. Cartouche automatique de contrôles de gaz du sang/CO-oxymètre, chimie Stat Profile Prime Plus® VET avec créatine, Stat Profile Prime Plus® VET Auto-Kassette mit Kreatinin für Blutgas-, CO-Oximeter- und Blutchemiekontrollen, Αυτόματο φύσιγγιο με κρεατινίνη επιπέδου ελέγχου χημείας, αερίων αίματος, CO-Οξυμέτρο Stat Profile Prime Plus® VET. Cartușcia con creatinina per controlli automatici chimici per gas ematici/CO-ossimetria Stat Profile Prime Plus® VET. Cartucho automático de controles de química, de CO-oxímetro e de gás no sangue Stat Profile Prime Plus® VET com creatinina, Stat Profile Prime Plus® VET vérégás, CO-oximéter, kémiai kontrollok automatikus patron kreatininnel, (מוסגרת בקרת כימיה אוטומטית עם קרטינין), (Stat Profile Prime Plus® VET 혈액 가스, CO-oks 측정, 화학 조절제 자동 카트리지, Stat Profile Prime Plus® VET 血気、一酸化炭素血気、化学制御液自動試劑盒 (含肌酐))

LOT 24207065 2025-11-29

CONTROL 1 2 3 4 5

Table with columns for Analyte, Control 1, Control 2, Control 3, Control 4, and Control 5. Rows include pH, H+, PCO2, PCO2, PO2, PO2, SO2, Hct, Na+, K+, Cl-, iCa, iCa, iMg, iMg, Glu, Glu, Lac, Lac, BUN, BUN, Urea, Urea, Creatinine, Creatinine, Creatinine, HbF, Hb, Hb, Hb, Hb, O2Hb, COHb, MetHb, HHb, tBil, tBil, tBil.

EN

Product Description: Aqueous quality control material for monitoring the performance of pH, PCO2, PO2, SO2, hematoct (Hct), fetal hemoglobin (HbF), total hemoglobin (HbT), oxymethemoglobin (COHb), methemoglobin (MetHb), and deoxyhemoglobin (iHb) in Levels 1, 2 and 3 as well as Na+, K+, Cl-, Ca, Mg, Glucose, Lactate, BUN (urea), and Creatinine in Level 4 and 5. A charge kit consisting of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation.

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Warnings and Cautions: DO NOT FREEZE. Mix the cartridge by gently inverting for several seconds. DO NOT SHAKE CARTRIDGE. Refer to Stat Profile Prime Plus VET Analyzer Instructions for Use Manual for complete information.

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Expected Ranges: The expected range for each parameter was determined at Nova Biomedical using replicate determinations on Nova analyzers. The expected range indicates the maximum deviations from the Mean Value that may be expected under differing laboratory conditions for instruments operating within specifications. Refer to Expected Ranges Table.

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200 Prospect Street Waltham, MA 02454 U.S.A.