Stat Profile® pHox®/Basic Controls Auto-Cartridge

Kontroll-Auto-Kassetten, Αυτόματη Κασέτα Διαλυμάτων Ελέγχου, Cartucho automático para controles, Cartouche CQ, Cartuccia automatica controlli, Cartucho automático de control, Kontroller, autokassett

CONTROL 1 2 3

LOT 24211021 2025-12-13

Table with columns for Analyte, Unit, and three Control Ranges (CONTROL 1, CONTROL 2, CONTROL 3). Rows include pH, H+, pCO2, pCO2, pO2, pO2, SO2, Hct, Hb, and Hb.

Product Description: Consists of 3 flexible bags within a cardboard carton. Each bag contains an aqueous quality control material for monitoring the measurement of pH, pCO2, pO2, SO2, Hematocrit (Hct), and Hemoglobin (Hb).

Controls: 1) Acidosis, with Low SO2, Low-Normal Hct/Hb; 2) Normal pH; 3) Alkalosis, with SO2, High Hct/Hb.

Intended Use: For in vitro diagnostic use for monitoring the performance of Nova Biomedical Stat Profile pHox/Basic Analyzers.

Methodology: Refer to Stat Profile pHox/Basic Analyzer Instructions For Use Manual.

Composition: A buffered bicarbonate solution, each control with a known pH. Solutions are equilibrated with known levels of O2, CO2, and N2. The reflectance characteristics give a signal that is equivalent to a known oxygen saturation value in whole blood.

Warnings and Cautions: Must be stored at 24-26°C for at least 24 hours before use. Intended for in vitro diagnostic use. Refer to the Stat Profile pHox/Basic Analyzer Instructions for Use for complete directions for use, including recommendations for use of controls, troubleshooting information, and methodology and Principles of the testing procedures.

Storage: Store at 2-8°C; DO NOT FREEZE. Expiration date is printed on each cartridge.

Directions for use: Cartridges must be stored at approximately 24-26°C for at least 24 hours prior to opening. Refer to Analyzer Instructions for Use for complete instructions. Verify that the lot number appearing on the Expected Ranges Table is the same lot number indicated on the control cartridge.

Limitations: PO2 values vary inversely with temperature (approximately 1%/°C). The Expected Range values are specific for instruments and calibrators manufactured by Nova Biomedical.

Traceability of Standards: Analytes are traced to NIST Standard Reference Materials.

Reference Intervals: Concentrations are formulated to represent three pH levels (Acidosis, Normal pH, and Alkalosis).

Users may wish to determine MEAN VALUES and EXPECTED RANGES in their own laboratory.*

*How to Define and Determine Reference Intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13

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EC REP Nova Biomedical GmbH Henssenring 13A, Geb. G 64546 Mörfelden-Walldorf Germany
Waltham, MA 02454 U.S.A.