Stat Profile Prime Plus® Blood Gas, CO-Oximeter, Chemistry Controls Auto-Cartridge with Creatinine

Cartucho automático con creatinina para controles de química, gases en sangre y cooxímetro Stat Profile Prime Plus®, Cartouche automatique de contrôles de gaz du sang/CO-oxymètre, chimie Stat Profile Prime Plus® avec créatine, Stat Profile Prime Plus® Auto-Kassette mit Kreatinin für Blutgas-, CO-Oximeter- und Blutchemiekontrollen, Αυτόματo φυσικόγιο υλικού ελέγχου χημείας ερωπιών αίματος, CO-Οξυμέτρο Stat Profile Prime Plus®, Cartuccia con creatinina per controlli automatici chimici per gas ematici/CO-ossimetria Stat Profile Prime Plus®, Cartucho automático de controles de química, de CO-oxímetro e de Gás no sangue Stat Profile Prime Plus® com creatinina, Stat Profile Prime Plus® vérgáz, CO-oximéter, kémiai kontrollok automatikus patron kreatininnel, מדידת כימיה אוטומטית עם קריאטינין Stat Profile Prime Plus® Co-Oximeter בדידת כימיה אוטומטית עם קריאטינין Stat Profile Prime Plus®, 血液ガス, CO オキシメーター, 生化学検査用コンロール自動カートリッジ(クレアチニン), 크레아티닌 사용 Stat Profile Prime Plus® 혈액 가스, CO-산소 농도계, 화학 조절자 자동 카트리지, Stat Profile Prime Plus® 血気 - 一酸化炭素血気儀 - 化学对照溶液自動试剂盒 (含肌酐)

LOT 24207065

2025-11-29

CONTROL 1 2 3 4 5

Table with columns for Expected Ranges, Rangos esperados, Plages attendues, Erwartungsbereiche, Αναμενόμενα εύρη, Intervalli previsti, Intervalos previstos, Várt tartományok, הצפויים הטווחים, 予測範囲, 예상 범위, 预期范围. Rows include pH, H+, PCO2, PO2, SO2, Hct, Na+, K+, Cl-, iCa, iMg, iGlu, Lac, BUN, Urea, Creatinine, HbF, tHb, O2Hb, MetHb, tBil, etc.

EN Product Description: Agrees quality control material for monitoring the performance of pH, PCO2, PO2, SO2, hematoct (Hct), fetal hemoglobin (HbF), total hemoglobin (Hb), total bilirubin (tBil), cyanmethemoglobin (COHb), methemoglobin (MetHb), and desoxyhemoglobin (HbD) in levels 1, 2 and 3 as well as Na+, K+, Cl-, Ca, Mg, Glucose, Lactate, BUN (urea), and Creatinine in Level 4 and 5. A charge kit consisting of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation. For use with Stat Profile Prime Plus Analyzers ONLY.

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Limitations: PO2 values vary inversely with temperature (approximately 1%/°C). Therefore, it is critical to follow the temperature guidelines described in "Directions for Use." The Expected Range values are specific for instruments and controls manufactured by Nova Biomedical. Once installed, each Stat Profile Prime Plus Cartridge may be used for a maximum of 21 days from the initial installation date on the system at which time the system will indicate the cartridge is invalid. Each cartridge may be inserted and removed from the analyzer a maximum of 6 times.