REF 48559, 41549

IVD 15°C 30°C CE

Stat Profile® pHox Ultra/Critical Care Xpress CO-Oximeter Calibrator Cartridge with Bilirubin and Deproteinizing Solution

Kalibrator-Patrone mit Bilirubin und Deproteinisierungslösung, Πακέτο βαθμονομητή με χοληρυθρίνη και διάλυμα αποπρωτεϊνισμού, Cartucho de calibrador con bilirubina y solución desproteinizante, Cartouche de calibrateur avec bilirubine et solution de déprotéinisation, Cartuccia calibratore con bilirubina e soluzione desproteinizzante, Cartucho de calibrador com solução de desproteinização e bilirubina, ビルリブリンおよび除蛋白液付キャリブレーション・カートリッジ, Calibrator Cartridge와 빌리루빈 및 단백질 제거액, 含有胆红素和去蛋白溶液的定标品盒

LOT 24180052 2025-06-28

Table with 4 columns: CALIBRATOR (≥70mL), DIL (≥180mL), SOLN F (≥620mL), and Deproteinizing Solution (≥90mL). Includes units (g/dL, g/L, mmol/L) and a WASTE CONTAINER warning.

Description: For use in calibrating the CO-Oximeter on the Stat Profile pHox Ultra/Critical Care Xpress (COX) Analyzer. The above table lists the concentrations of the constituents for the calibrator, flush, deproteinizing and diluent solutions contained in the Calibrator Cartridge.

Intended Use: For the quantitative determination of total hemoglobin, oxyhemoglobin, carboxyhemoglobin, methemoglobin, deoxyhemoglobin and total bilirubin in human blood using the Stat Profile pHox Ultra/Critical Care Xpress CO-Oximeter Analyzer.

Methodology: Refer to the Stat Profile pHox Ultra/Critical Care Xpress Analyzer Instructions for Use Manual for Test Methodologies and Principles.

Composition: Solution volumes are listed in the table above. Contains no constituents of human origin, however good laboratory practice should be followed during handling of these materials.

Warnings and Cautions: Intended for in vitro diagnostic use. Refer to the Stat Profile pHox Ultra/Critical Care Xpress Analyzer Instructions for Use Manual for complete directions. Follow standard practices for handling laboratory reagents.

Storage: Store at 15°C - 30°C. Expiration date is printed on the front of the cartridge.

Directions for Use: Refer to the Stat Profile pHox Ultra/Critical Care Xpress Analyzer Instructions for Use Manual for complete directions for installation of the cartridge. Cartridge is ready for use.

Limitations: For use on Nova Biomedical Analyzers only. If not stored at room temperature and out of direct sunlight, the solution may deteriorate.

Traceability of Standards: Traceable to precise hemoglobin measurement using the reference cyanmethemoglobin method.

Reference: ICSH, Recommendations for Haemoglobinometry in Human Blood. British Journal of Haematology, 13 (suppl.) pp 71-75, 1967.

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EC REP Nova Biomedical GmbH Hessenring 13A, Geb. G 64546 Mörfelden-Walldorf Germany
200 Prospect Street Waltham, MA 02454 U.S.A. www.novabiomedical.com