LOT 24221050

2026-08-02

CONTROL 1 2 3

IVD 15°C 30°C

Table with 5 columns: Parameter, Unit, Control 1 (min - x̄ - max), Control 2 (min - x̄ - max), Control 3 (min - x̄ - max). Rows include pH, H+, PCO2, PO2, SO2, Hct, Hb, and Hb.

Product Description

An aqueous quality control material for monitoring the measurement of pH, pCO2, pO2, SO2, hematoct (Hct), hemoglobin (Hb) for use with Nova Biomedical analyzers ONLY. Formulated at three levels:

- CONTROL 1 Acidosis, Low SO2, Low-Normal Hct/Hb
CONTROL 2 Normal
CONTROL 3 Alkalosis, with SO2, Normal-High Hct/Hb

Intended Use

For in vitro diagnostic use for monitoring the performance of Stat Profile pHox/Basic Analyzers.

Methodology

Refer to Stat Profile pHox/Basic Analyzer Instructions For Use Manual.

Composition

A buffered bicarbonate solution, each control with a known pH. Solutions are equilibrated with known levels of O2, CO2, and N2. The reflectance characteristics give a signal that is equivalent to a known oxygen saturation value in whole blood.

Warnings and Cautions:

Must be stored at 24-26°C for at least 24 hours before opening. Intended for in vitro diagnostic use. Refer to the Stat Profile pHox/Basic Analyzer Instructions for Use Manual for complete directions for use.

Storage

Store at 15-30°C. DO NOT FREEZE. Expiration date is printed on each ampule container. Do not store unused opened ampules.

Directions for use

Contents must be shaken well prior to use. Controls must be stored at approximately 24-26°C for at least 24 hours prior to opening. Refer to Analyzer Instructions for Use Manual for complete instructions.

Limitations

PO2 values vary inversely with temperature (approximately 1%/°C).

The Expected Range values are specific for analyzers and calibrators manufactured by Nova Biomedical.

Traceability of Standards

Analyzers are traced to NIST Standard Reference Materials.

Reference Intervals

Concentrations are formulated to represent three pH levels (Acidosis, Normal pH, and Alkalosis).

The expected clinical range of values for these analytes in patient blood is referenced in Tietz, NW ed 1986 Textbook of Clinical Chemistry, W.B. Saunders Co.

Users may wish to determine MEAN VALUES and EXPECTED RANGES in their own laboratory!

Expected Ranges

The EXPECTED RANGE for each analyte was determined at Nova Biomedical by using multiple runs of each level of control at 37°C on multiple analyzers. The EXPECTED RANGE indicates the maximum deviations from the mean value that may be expected under differing laboratory conditions for analyzers operating within specifications.

\*How to Define and Determine Reference Intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13

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Waltham, MA 02454 U.S.A. www.novabiomedical.com