Linearity Standard Set A

Lineariitätsstandard Satz A, Πρότυπο Γραμμικότητας Set A, Estándar de linealidad Set A, Set A standard de linearité, Standard di linearità Set A, Conjunto A de padrão de linearidade, Linearity Standard Set A

LEVEL 6

Expected Ranges, Erwartete Bereiche, Αναμενόμενα Εύρη, Rangos esperados, Fourchettes attendues, Intervalli previsti, Gamas previstas, Förväntade områden			
LOT 24208023		2026-01-26	LEVEL 6
NOVA 1-16	Na ⁺	mmol/L	90 - 110
	K ⁺	mmol/L	90 - 110
	Cl ⁻	mmol/L	180 - 220

Product Description
Aqueous solution containing Na⁺, K⁺, and Cl⁻. For use on NOVA analyzers ONLY.

Intended Use
For in vitro diagnostic use with NOVA Analyzers to verify calibration, analytical linearity, estimate test imprecision, and detect systematic analytical deviations that may arise from calibrator cartridge or analytical instrument variation.

Methodology
Refer to Nova Analyzer Instructions for Use for Methodology and Principles of the testing procedures.

Composition
Buffered solutions containing Na⁺, K⁺, and Cl⁻. Minimum volume 50 mL. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF. NCCLS DOCUMENT M29-T2).

Warnings and Cautions
Intended for In Vitro Diagnostic Use. Clumping or particulate matter in the solution is an indication of possible deterioration. DO NOT FREEZE. Follow standard practices for handling laboratory reagents. Discard in accordance with local requirements.

Storage
Store at 15-30°C. After opening, the solution is stable for 6 months.

Directions for Use
Linearity Standards are recommended for use as frequently as required by local regulatory and hospital requirements.

Limitations
It is recommended that each laboratory establish the upper and lower limits of its reportable range for each sample type and devise a policy for processing specimens that fall outside of this range.

Traceability of Standards
Analytes traced to NIST Standard Reference Materials.

Expected Ranges
The expected range indicates the maximum deviations from the mean value that may be expected under differing laboratory conditions for instruments operating within specifications. The expected ranges for the analytes listed in the table above were verified using replicate determinations on Nova analyzers.

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