
Stat Profile® Prime VET CCS Comp Calibrator Cartridge

Kalibratorkassette, Καρέττα Βαθμονομητή, Cartucho del calibrador, Cartouche étalon, Cartuccia calibratore, Calibrator Cartridge, Kalibratorokassett, Kalibrátorkazetta, 고정기 카트리지, キャリブレータカートリッジ, 校准液盒

LOT 24228026

2026-02-06

| CAL A ≥ 500 mL | | CAL B ≥ 250 mL | | CAL F ≥ 720 mL | | SOLN R ≥ 300 mL | |
|------------------|--------------|-----------------|-------------|------------------|-----------|---|---------|
| pH | 7.359 | pH | 6.840 | PCO ₂ | 51.6 mmHg | KCl | 2 mol/L |
| PCO ₂ | 24.4 mmHg | PO ₂ | 167.0 mmHg | PO ₂ | 99.2 mmHg |  | |
| Na ⁺ | 133.0 mmol/L | Na ⁺ | 68.5 mmol/L | | | | |
| K ⁺ | 4.0 mmol/L | K ⁺ | 10.0 mmol/L | | | | |
| Cl ⁻ | 106.0 mmol/L | Cl ⁻ | 46.0 mmol/L | | | | |
| iCa | 1.10 mmol/L | iCa | 2.20 mmol/L | | | | |
| Glu | 80 mg/dL | Glu | 200 mg/dL | | | | |
| Lac | 2.0 mmol/L | Lac | 10.0 mmol/L | | | | |

EN
Product Description
 The above table lists the concentrations of constituents for the calibrators and solutions contained in the Calibrator Cartridge. Contains the following aqueous solutions within individual flexible bags housed in a cardboard box: 3 Calibrators (Calibrator A, B, and F), reference solution, and a flexible waste bag. Each bag includes a fihnent with a septa. The exposed bag fihnents are arranged in a line along the rear of the cartridge box. These septa are pierced during the insertion of the cartridge into the analyzer.

Intended Use
 For the quantitative determination of Na, K, Cl, iCa, Glucose, Lactate, pH, PCO₂, PO₂, and Hct in blood using the Stat Profile Prime VET CCS Comp Analyzer.

Methodology
 Refer to Prime Vet CCS Comp Analyzer Instructions For Use Manual for Test Methodologies and Principles.

Composition
 Calibrator volumes are listed in the table above. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF: NCCLS DOCUMENT M29-T2).

Warnings and Cautions:
 For Veterinary Use Only. Not for Human Use. Refer to the Prime Vet CCS Comp Analyzer Instructions for Use Manual for complete directions for use, installation instructions, recommendations for calibration intervals, Troubleshooting information, Methodology and reagents. Discard Calibrator Cartridges and Waste products in accordance with local guide lines. Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

Storage
 Store at 15-30°C. Expiration date is printed on the front of the Calibrator Cartridge.

Directions For Use
 Refer to the Prime Vet CCS Comp Analyzer Instructions for Use for complete directions for installation of the cartridge, and recommendations for Calibration Intervals. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents. Cartridge is ready for use.

Traceability of Standards
 Analytes traced to NIST Standard Reference Materials.

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