**Stat Profile Prime Plus® Blood Gas, CO-Oximeter, Chemistry Controls Auto-Cartridge with Creatinine**

Cartucho automático con creatinina para controles de química, gases en sangre y cooxímetro Stat Profile Prime Plus®, Cartouche automatique de contrôles de gaz du sang/CO-oxymètre, chimie

Stat Profile Prime Plus® avec créatine, Stat Profile Prime Plus® Auto-Kassette mit Kreatinin für Blutgas-, CO-Oximeter- und Blutchemiekontrollen, Αυτόματο φυσικό υλικό ελέγχου χρησίας αερίων αίματος, CO-Οξύμετρο Stat Profile Prime Plus®, Cartuccia con creatinina per controlli automatici chimici per gas ematici/CO-ossimetro Stat Profile Prime Plus®, Cartucho automático de controles de química, de CO-oxímetro e de Gás no sangue Stat Profile Prime Plus® com creatinina, Stat Profile Prime Plus® vérágz, CO-oximéter, kémiai kontrollok automatikus patron

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FR

## Product Description

Acceso control material for monitoring the performance of pH, PCO<sub>2</sub>, PO<sub>2</sub>, SO<sub>2</sub>, hematocrit (Hct), fetal hemoglobin (HbF), total hemoglobin (Hb), total bilirubin (TBIL), oxyhemoglobin (Hb), carboxyhemoglobin (COHb), methemoglobin (MetHb) and deoxyhemoglobin (DeoxyHb) in Levels 1, 2 and 3 as well as Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, Ca<sup>2+</sup>, Mg<sup>2+</sup>, Glucose, Lactate, BUN (Urées), and Creatinine in Levels 4 & 5. The cartridge consisting of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation. For use with Stat Profile Prime Plus Analyzers ONLY.

## Intended Use

Designed for *in vitro* diagnostic use by healthcare professionals for monitoring the performance of the Stat Profile Prime Plus Analyzers.

## Methodology

Refer to Stat Profile Prime Plus Analyzer Instructions for Use for Methodology and Principles.

## Composito

Composito Levels 1, 2 and 3 are buffered bicarbonate solutions containing dye, salts and preservatives. Each level has a known pH and is supplemented with a known O<sub>2</sub>, CO<sub>2</sub> and Na<sup>+</sup>. Control Level 4 and 5 are buffered solutions containing known concentrations of Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, Ca<sup>2+</sup>, Mg<sup>2+</sup>, Glucose, Lactate, BUN (Urées) and Creatinine.Components Levels 1, 2 and 3 are buffered bicarbonate solutions containing dye, salts and preservatives. Each pouch contains a minimum of 100 mL. Components Levels 4 and 5 are buffered solutions containing known concentrations of Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, Ca<sup>2+</sup>, Mg<sup>2+</sup>, Glucose, Lactate, BUN (Urées) and Creatinine.

## Warnings and Cautions

NO FREEZING. Mix the cartridge by gently inverting for several seconds. DO NOT SHAKE CARTRIDGE. Refer to Stat Profile Prime Plus Analyzer Instructions for Use for complete information.

Intended for *in vitro* diagnostic use. Follow standard practices for handling laboratory reagents.

## Storage

Store at 2-4°C (37-46°F). DO NOT FREEZE.

## Directions for use

Ensure controls are room temperature prior to installation.

Level 4 and Level 5 must be charged with creatinine syringes prior to installation of the Control Cartridge. The 2 syringes are labeled and attached to the syringes on the cartridge.

Activate the cartridge as follows:

1. Hold the syringe with tip side down and gently invert.

2. Allow the air bubbles to rise to the top of the syringe. Remove the protective cover from the needle.

3. Match the color and lase of the syringe to the appropriate lase and insert needle.

4. Slowly depress syringe plunger until the contents are dispensed. DO NOT PULL BACK ON THE PLUNGER TO FLUSH CARTRIDGE.

5. Remove needle/syringe assembly from lase and discard in an appropriate sharps container.

6. Repeat Steps 1-4 for the next Control.

7. Mix the cartridge by gently inverting for 1 minute. Cartridge is ready for use.

Verify that the Lot Number on the Expected Range Table corresponds to the Lot Number on the cartridge. Refer to Stat Profile Prime Plus Analyzer Instructions for Use for complete directions.

## Limitations

PO<sub>2</sub> values vary inversely with temperature (approximately -0.1%/°C). Therefore, it is critical to follow the temperature guidelines described in "Directions for use." The Expected ranges are specific for instruments and controls manufactured by Nova Biomedical. Carboxyhemoglobin (COHb) and Methemoglobin (MetHb) are traceable using Spectrophotometry. Analyses are traced to NIST Standard Reference Materials.

## Reference Intervals

The expected range for each parameter was determined at Nova Biomedical using replicate determinations on Nova analyzers. The expected range indicates the maximum deviations from the Mean Value that may be expected under differing laboratory conditions for instruments operating with specifications. Refer to Expected Ranges table.

## Traceability of Standards

Total hemoglobin (Hb) and Methemoglobin (MetHb) are traceable using Cyanmethemoglobin method. Carboxyhemoglobin (COHb) and Oxyhemoglobin (O<sub>2</sub>Hb) are traceable using Spectrophotometry. Analyses are traced to NIST Standard Reference Materials.

## Expected Ranges

The expected range for each parameter was determined at Nova Biomedical using replicate determinations on Nova analyzers. The expected range indicates the maximum deviations from the Mean Value that may be expected under differing laboratory conditions for instruments operating with specifications. Refer to Expected Ranges table.

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