REF 48903, 48904

IVD 20°C - 80°C CE

Stat Profile® pHox Ultra/Critical Care Xpress Blood Gas Controls Auto-Cartridge

Blutgas-Kontroll-Auto-Kassetten, Αυτόματη Κασέτα Διαλυμάτων Ελέγχου Αερίων Αίματος, Cartucho automático para controles de gases en sangre, Cartouche CQ pour les gaz du sang, Cartuccia automatica controlli, Cartucho automático de controlo de gases no sangue, Blodgaskontroller, autokasset

CONTROL 1 2 3

LOT 24226037 2026-01-25

	Expected Ranges, Erwartete Bereiche, Αναμενόμενα Εύρη Rangos esperados, Fourchettes attendues, Intervalli previsti, Gamas previstas, Förväntade områden		
	CONTROL 1	CONTROL 2	CONTROL 3
	min - \bar{x} - max	min - \bar{x} - max	min - \bar{x} - max
pH	7.159 - 7.184 - 7.209 ✓	7.366 - 7.391 - 7.416 ✓	7.518 - 7.545 - 7.572 ✓
H ⁺	nmol/L 69.41 - 65.53 - 61.86 ✓	43.03 - 40.62 - 38.35 ✓	30.32 - 28.49 - 26.77 ✓
pCO ₂	mmHg 56.1 - 61.4 - 66.6 ✓	40.0 - 43.5 - 47.0 ✓	21.5 - 24.0 - 27.0 ✓
pCO ₂	kPa 7.45 - 8.16 - 8.85 ✓	5.32 - 5.78 - 6.25 ✓	2.87 - 3.20 - 3.60 ✓
pO ₂	mmHg 50.1 - 58.1 - 66.1 ✓	99.7 - 107.7 - 115.7 ✓	125.6 - 137.6 - 149.6 ✓
pO ₂	kPa 6.66 - 7.73 - 8.79 ✓	13.26 - 14.32 - 15.39 ✓	16.70 - 18.30 - 19.89 ✓
SO ₂	% 54.1 - 59.1 - 63.1 ✓		88.7 - 92.7 - 96.7 ✓
Hct	% 32 - 34 - 37 ✓		46 - 49 - 52 ✓
Hb	g/dL 9.6 - 11.1 - 12.6 ✓		14.8 - 16.3 - 17.8 ✓
Hb	mmol/L 5.96 - 6.89 - 7.82 ✓		9.19 - 10.12 - 11.05 ✓

Product Description
An aqueous quality control material for monitoring the measurement of pH, pCO₂, pO₂, SO₂, hematox (Hct), and hemoglobin (Hb) for use with Nova Biomedical analyzers ONLY. Each controls Auto-Cartridge contains controls formulated at three clinically significant levels:

Normal pH
Acidosis, with low SO₂, low-normal Hct/Hb
Alkalosis, with SO₂, High Hct/Hb

Intended Use
For in vitro diagnostic use for monitoring the performance of Nova Biomedical Stat Profile pHox Ultra/Critical Care Xpress Analyzers.

Methodology
Refer to Stat Profile pHox Ultra/CCX Analyzer Instructions For Use Manual.

Composition
A buffered bicarbonate solution, the Controls Auto-Cartridge consists of 3 flexible bags, each with a known pH. Solutions are equilibrated with known levels of O₂, CO₂, and N₂. The reflectance characteristics give a signal that is equivalent to a known oxygen saturation value in whole blood. The conductivity characteristics are equivalent to a known hematocrit value in whole blood. The reflectance and conductivity signal is equivalent to a known hemoglobin value in blood. Mold inhibited. Each control flexible pouch contains a minimum value of 100mls. Contains no constituents of human origin, however good laboratory practice should be followed during handling of these materials. (REF: NCCLS DOCUMENT M23-T2)

Warnings and Cautions:
Must be stored at 24-26°C for at least 24 hours before using. DO NOT SHAKE the cartridge. Intended for in vitro diagnostic use. Refer to the Stat Profile pHox Ultra/Critical Care Xpress Analyzer Instructions for Use for complete directions for use, recommendations for use of controls, troubleshooting information, and methodology and Principles of the testing procedures. Follow standard practices required for handling laboratory reagents.

Storage
Store at 2-8°C. DO NOT FREEZE. It is critical to follow the room temperature equilibrium instructions prior to use as described in "Directions for Use". Each cartridge has a lot number and expiration date printed on the label.

Directions for use
Controls must be stored at approximately 24-26°C for at least 24 hours prior to opening. Refer to Analyser Instructions for Use for complete instructions. Verify that the lot number appearing on the Expected Ranges Table is the same lot number indicated on the control material being analyzed. The Auto-Cartridge should be mixed by gently inverting the pack for several seconds. DO NOT SHAKE the cartridge.

Limitations
PO₂ values vary inversely with temperature (approximately 1%/°C).
The Expected Range values are specific for instruments and controls manufactured by Nova Biomedical. Once installed, each Stat Profile pHox Ultra/CCX Auto-Cartridge may be used for a maximum of 35 days from the initial installation date on the system at which time the system will indicate the cartridge is invalid.

Traceability of Standards
Analytes are traced to NIST Standard Reference Materials.

Reference Intervals
Concentrations are formulated to represent three pH levels (Acidosis, Normal pH, and Alkalosis), and with Low-Normal and High Hct/Hb levels.

The expected clinical range of values for these analytes in patient blood is referenced in Tietz, NW ed 1986 Textbook of Clinical Chemistry, W.B. Saunders Co.

Users may wish to determine MEAN VALUES and EXPECTED RANGES in their own laboratory!

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