
REF 53104, 53360, 53359, 53463, 53464, 53465

# Stat Profile® Prime Calibrator Cartridge ABG

Kalibratorkassette, Κάσета Βαθμονομητή, Cartucho del calibrador, Cartouche étalon, Cartuccia calibratore, Calibrator Cartridge, Kalibratorkassett

LOT 24247063

2026-02-12

CAL A	≥ 500 mL	CAL B	≥ 250 mL	CAL F	≥ 720 mL	SOLN R	≥ 300 mL
pH	7.351	pH	6.840	PCO <sub>2</sub>	50.7 mmHg	KCl	2 mol/L
PCO <sub>2</sub>	24.7 mmHg	PO <sub>2</sub>	167.0 mmHg	PO <sub>2</sub>	99.0 mmHg	<b>WASTE CONTAINER</b> W	

**Product Description**  
The above table lists the concentrations of constituents for the calibrators and solutions contained in the Calibrator Cartridge. Contains the following aqueous solutions within individual flexible bags housed in a cardboard box: 3 Calibrators (Calibrator A, B, and F), reference solution, and a flexible waste bag. Each bag includes a fihnent with a septa. The exposed bag fittings are arranged in a line along the cartridge box. These septa are pierced during the insertion of the cartridge into the analyzer.

**Intended Use**  
For the quantitative determination of pH, PCO<sub>2</sub>, and PO<sub>2</sub> in human blood using the Nova Biomedical Prime ABG Analyzer.

**Methodology**  
Refer to Prime Analyzer Instructions For Use Manual for Test Methodologies and Principles.

**Composition**  
Calibrator volumes are listed in the table above. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF. NCCLS DOCUMENT M29-T2).

**Warnings and Cautions:**  
Intended for *in vitro* diagnostic use.  
Refer to the Prime Analyzer Instructions for Use Manual for complete directions for use, installation instructions, recommendations for calibration intervals, Troubleshooting information, Methodology and Principles of the testing procedure. Follow standard practices required for handling laboratory reagents. Discard Calibrator Cartridges and Waste products in accordance with local guidelines. Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

**Storage**  
Store at 15-30°C. Expiration date is printed on the front of the Calibrator Cartridge. Directions For Use Refer to the Prime Analyzer Instructions for Use for complete directions for installation of the cartridge, and recommendations for Calibration Intervals. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents. Cartridge is ready for use.

**Traceability of Standards**  
Analytes traced to NIST Standard Reference Materials.

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