REF 49187

IVD 15°C 30°C

Stat Profile® pHox Ultra/Critical Care Xpress CO-Oximeter Controls

CO-Oximeter-Kontrollen, Διαλύματα Ελέγχου Οξυμέτρου-CO, Controles de CO-Oxímetro, Contrôles par oxymètre CO, Controlli qualità co-ossimetro, Controles de co-oxímetro, CO-oximeter-kontroller, CO-酸素濃度計コントロール, CO-산소 농도계 조절, CO-血氧計质控品

CONTROL 7 8 9

LOT 24239029

2026-08-15

LOT

Expected Ranges, Erwartete Bereiche, Αναμενόμενα εύρη, Rangos esperados, Fourchettes attendues, Intervalli previsti, Gamas previstas, Förväntade områden, 予測範囲, 예상 범위, 预期范围

		CONTROL 7	CONTROL 8	CONTROL 9
		min - \bar{x} - max	min - \bar{x} - max	min - \bar{x} - max
tHb	g/dL	6.0 - 6.7 - 7.4	15.1 - 16.1 - 17.1	20.5 - 22.0 - 23.5
tHb	g/L	60.0 - 67.0 - 74.0	151.0 - 161.0 - 171.0	205.0 - 220.0 - 235.0
tHb	mmol/L	3.7 - 4.2 - 4.6	9.4 - 10.0 - 10.6	12.7 - 13.6 - 14.6
O ₂ Hb	%	78.7 - 83.7 - 88.7	47.1 - 51.1 - 55.1	20.9 - 23.4 - 25.9
COHb	%	4.2 - 8.2 - 12.2	15.8 - 19.8 - 23.8	23.9 - 27.9 - 31.9
MetHb	%	3.6 - 6.6 - 9.6	15.5 - 18.5 - 21.5	25.9 - 28.9 - 31.9
HHb	%	0.0 - 4.0 - 8.0	6.6 - 10.6 - 14.6	15.8 - 19.8 - 23.8
tBil	mg/dL	3.6 - 4.0 - 4.4	8.2 - 10.2 - 12.2	16.7 - 20.7 - 24.7

Description
An aqueous quality control material for monitoring the measurement of total hemoglobin (tHb), oxyhemoglobin (O₂Hb), carboxyhemoglobin (COHb), methemoglobin (MetHb), deoxyhemoglobin (HHb) and total bilirubin (tBil). For use with Stat Profile pHox Ultra/Critical Care Xpress (CCX) Analyzers ONLY. The controls are formulated at three levels.

Intended Use
For in vitro diagnostic use by health care professionals for monitoring the performance of Stat Profile pHox Ultra/CCX Analyzers.

Methodology
Refer to Stat Profile pHox Ultra/CCX Analyzer's Instructions for Use Manual for use, Methodology and Principles.

Composition
Contains dyes, viscosity enhancers and a preservative in an aqueous salt solution. Each ampule contains 1.7 mL. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials (REF. NCCLS DOCUMENT M29-T2).

Warnings and Cautions
Intended for in vitro diagnostic use. Refer to the Stat Profile pHox Ultra/CCX Analyzer's Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory reagents. Once ampule is opened, discard unused portion in accordance with local guidelines.

Storage
Store at 15-30°C. Exposure to direct sunlight should be avoided. Expiration date is printed on each ampule. Do not store unused opened ampules.

Directions for use
Refer to the analyzer's Instructions for Use Manual for complete directions. Each ampule is intended for single use only. However, multiple analyses can be processed from one ampule. Verify that the lot number appearing on the Expected Ranges Table is the same lot number indicated on the control ampule.

Limitations
For use on the Stat Profile pHox Ultra/CCX Analyzers only. If not stored at room temperature and out of direct sunlight, the solution may deteriorate. Indications of deterioration may be high or low range control values or a noticeable change in calibration or both. Cloudiness or particulate matter may be a sign of possible deterioration.

This product is warranted to perform as described in this product insert, and Nova Biomedical disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall Nova Biomedical be liable for any consequential damages arising out of the aforesaid expressed warranty.

Traceability of Standards
Total Hemoglobin (tHb) and Methemoglobin (MetHb) are traceable by using Cyanmethemoglobin method. Carboxyhemoglobin (COHb) and Oxyhemoglobin (O₂Hb) are traceable using Spectrophotometry.

Reference Intervals
Concentrations are formulated to represent three levels.

Users may wish to determine MEAN VALUES and EXPECTED RANGES in their own laboratory.¹

Expected Ranges
The EXPECTED RANGE for each parameter was determined at Nova Biomedical by using multiple runs of each level of control on properly calibrated instruments. The EXPECTED RANGE indicates the maximum deviations from the mean value that may be expected under differing laboratory conditions for instruments operating within specifications. Refer to Expected Ranges Table.

¹How to Define and Determine Reference Intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13

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